

MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit

for use with the PrepSEQ™ Mycoplasma Sample Preparation Kit

Catalog Numbers 4460623, 4460626

Publication Number 4465874

Revision C.0



Manufacturer: Thermo Fisher Scientific | 7 Kingsland Grange | Warrington, Cheshire WA1 4SR | United Kingdom™

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4465874

Revision	Date	Description
C.0	24 May 2018	Updated template, legal, and content information. Reorganized content. Added information about using the AccuSEQ™ Software v2.0 Mycoplasma SEQ module.
B.0	04 Dec 2013	Baseline for this revision history.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

TRADEMARKS: All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2018 Thermo Fisher Scientific Inc. All rights reserved.

Contents

■	Product information	5
	Product description	5
	Contents and storage	5
	Required materials not supplied	6
	Workflow	7
■	Methods	8
	Prepare the sample	8
	Prepare for PCR using AccuSEQ™ Software v2.0 Mycoplasma SEQ module	8
	Create a new experiment	8
	Define the experiment properties	9
	Setup the samples and controls	9
	View the plate layout	10
	Save the experiment	11
	Prepare the kit reagents and premix solution	11
	Prepare the PCR reactions	12
	Start the run	13
	Analyze the results	13
	Set the baseline and threshold values	13
	Review the Results Summary	14
	Guidance for test samples	15
	Guidance for controls	15
	Guidance for inconclusive results with AccuSEQ™ software v2.0	15
	Example results with AccuSEQ™ Software v. 2.1.1	16
	Positive control	16
	Negative control	17
	Blank extraction control	18
	Positive extraction control	19
	Inhibition control and positive control	20
	Test sample: Negative result	21
	Test sample: Positive result	22
	Test sample: Positive result with decreased detection of DPC	23
	PCR inhibition	24
	Multicomponent plots	25

■ APPENDIX A	Use the kit with 7500 System SDS Software v1.4 or later	26
	Prepare the kit reagents and premix solution	26
	Prepare the PCR reactions	27
	Seal the plates	28
	Prepare the plate document	29
	Perform PCR	30
	Analyze the results	31
	Set the baseline and threshold values	31
	Guidance for test samples	32
	Guidance for controls	32
	Example positive results with SDS v1.4 software	33
	Example positive control extraction results with SDS v1.4 software	34
	Example negative results with SDS v1.4 software	35
■ APPENDIX B	Troubleshooting	36
	AccuSEQ™ 2.0 software	36
	MycoSEQ™ kit	37
■ APPENDIX C	Background information	39
■ APPENDIX D	Kit specificity	40
	Sensitivity	40
	Kit specificity	40
	Inclusivity – detectable species	40
	Exclusivity – undetectable organisms	41
■ APPENDIX E	Good PCR practices	43
	Good laboratory practices for PCR and RT-PCR	43
	Plate layout suggestions	43
	Documentation and support	44
	Related documentation	44
	Customer and technical support	45
	Limited product warranty	45



Product information





IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Product description

The MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit detects *Mycoplasma* species simply, reliably, and rapidly. To detect the presence of these microorganisms, the assay uses the polymerase chain reaction (PCR) to amplify a target unique to a wide variety of *Mycoplasma* species.

Contents and storage

Table 1 MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit (Cat. No. 4460623)

Contents ^[1]	Cap color	Amount ^[2]	Storage
Box 1: MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit			
10× <i>Mycoplasma</i> Real-Time PCR Primer Mix	 blue	325 µL	–25°C to –15°C on receipt 2–8°C after first use
Negative Control	 white	1,000 µL	
2× <i>PowerSYBR™</i> Green PCR Master Mix	 white	2 × 1,000 µL	–25°C to –15°C on receipt, protected from light. 2–8°C after first use, protected from light
Box 2: MycoSEQ™ Discriminatory Positive/Extraction Control			
MycoSEQ™ Discriminatory Positive/Extraction Control, 1,000 copies/µL	 yellow	700 µL	–25°C to –15°C

^[1] To purchase the MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit that includes the PrepSEQ™ Mycoplasma Sample Preparation Kit, use Catalog Number 4460626.

^[2] The kit contains reagents for 100 reactions.



Required materials not supplied


Unless otherwise indicated, all materials are available through **thermofisher.com**.
MLS: Fisher Scientific (**fisherscientific.com**) or other major laboratory supplier.

Item	Source
Instruments; choose one:	
7500 Fast Real-Time PCR System with AccuSEQ™ software v2.0 or later	Contact your local sales representative
[Optional] 7500 Real-Time PCR System	
Consumables	
Disposable gloves	MLS
Aerosol-resistant pipette tips	MLS
Pipettors: <ul style="list-style-type: none"> Positive-displacement Air-displacement Multichannel 	MLS
MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 20 plates, 0.2-mL well; for use with 7300, 7500, and 7900HT Fast Real-Time PCR Systems	4306737 ^[1]
MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL, 20 plates; for use with 7500 Fast Real-Time PCR System	4346906
MicroAmp™ Optical 96-Well Reaction Plate with Barcode & Optical Adhesive Films, 100 plates with covers; for use with 7300 and 7500 Fast Real-Time PCR Systems	4314320
MicroAmp™ Optical 8-Cap Strips, 300 strips	4323032
MicroAmp™ Optical Adhesive Film Kit, 20 covers, 1 applicator, 1 optical cover compression pad	4313663
MicroAmp™ Optical Adhesive Film, 25 or 100 covers	4360954, 25 covers 4311971, 100 covers

^[1] Not recommended for use with the 7500 Fast system. For 7500 Fast system reactions, use Cat. No. 4346906.



Workflow





Methods

IMPORTANT! This chapter describes how to prepare and run PCR samples using AccuSEQ™ Software v2.0. If you are using SDS software v1.4 or later, see Appendix A, “Use the kit with 7500 System SDS Software v1.4 or later”.

IMPORTANT! For information on how to avoid PCR contamination, see Appendix E, “Good PCR practices”.

Prepare the sample

Prepare the DNA template for the PCR reactions using the PrepSEQ™ *Mycoplasma* Nucleic Acid Extraction Kit.

For more information, see:

- The *PrepSEQ™ Sample Preparation Kits for Mycoplasma, MMV, and Vesivirus User Guide* (Pub. No. 4465957)
- The *PrepSEQ™ Express Nucleic Acid Extraction Kit for Mycoplasma, MMV, and Vesivirus Detection User Guide* (Pub. No. MAN0016799)

Prepare for PCR using AccuSEQ™ Software v2.0 Mycoplasma SEQ module

Create a new experiment

1. In the desktop, double-click the AccuSEQ™ software icon to start the software.
2. Log into the software. In the **Home** screen, click **Create MycoSEQ Experiment** to open the **Mycoplasma Assay v2.0** workflow.



Define the experiment properties

In the Experiment properties screen:

1. Enter an experiment name.

2. (Optional) Enter a barcode to identify the reaction plate.
3. (Optional) Enter comments to describe the experiment.
4. Verify the SEQ experiment type and assay to use, then click **Next**.

Setup the samples and controls

In the Sample Setup screen:

1. Specify the number of samples and replicates:

Sample Type	Sample Name	Name Fill	Plot Color	Number of Wells - Sample	Number of Wells - Inhibition Controls
UNKNOWN	Sample 1			1	1
UNKNOWN	Sample 2			1	1
UNKNOWN	Sample 3			1	1
UNKNOWN	Sample 4			1	1
UNKNOWN	Sample 5			1	1
Positive Control	POS 1			1	0
Negative Control	NEG 1			1	0



Field	Minimum entry ^[1]
Samples	1
Sample replicates	1
Inhibition control replicates for each sample	0
Positive control replicates	1
Negative control replicates	1

[1] We recommend that you use at least one negative and one positive control per run, and at least one inhibition control per sample.

2. Set the Sample volume to 10 µL per reaction.
3. Enter sample names, and *(optional)* set plot colors.
4. Click **Next**.

View the plate layout

The AccuSEQ™ software v2.0 uses the sample information that you enter in the **Sample Setup** screen to fill the wells in the plate layout and to calculate the required reaction component volumes for each sample type, based on the Mycoplasma Real-Time PCR Detection Kit guidelines.

To view and edit the plate layout before starting an instrument run:

1. Click **Setup ▶ Plate Layout** in the navigation pane.
2. Review the initial well selections in the Plate Layout screen. Drag-and-drop samples to create the layout of your choice.

Here is an example:

3. Review the **Sample Setup** window to ensure that the number of Unknowns, Inhibition Controls, Positive Controls, and Negative Controls match your experiment sample setup. In the example, this is 5 Unknowns, 5 Inhibition Controls, 1 Positive Control, and 1 Negative Control.
4. Review the run method and click **Next**.



Save the experiment

1. At the bottom of the AccuSEQ™ software screen, click **Save & Finish**.
2. In the **Save Experiment** dialog box, verify the Mycoplasma Presence Absence Detection by MycoSEQ Example.eds file name, then click **Save**.

Prepare the kit reagents and premix solution

1. Thaw all kit reagents completely.
2. Vortex briefly, then spin down the reagents.
3. Prepare the Premix Solution according to the following table.

Component for premix solution	Volume for one 30-μL reaction	Volume for four 30-μL reactions ^[1]
<i>Power</i> SYBR™ Green PCR Master Mix, 2×	15.0 μL	66.0 μL
<i>Mycoplasma</i> Real-Time PCR Primer Mix, 10×	3.0 μL	13.2 μL
Total premix solution volume	18.0 μL	79.2 μL

^[1] Includes 10% excess to compensate for pipetting errors.

4. Mix the Premix Solution by gently pipetting up and down, then cap the tube.



Prepare the PCR reactions

1. Dispense the following into each well to be used, gently pipetting at the bottom of the well.

To prepare...	In each tube or well...
Negative control reaction	<ul style="list-style-type: none">• Add 18 μL of Premix Solution• Add 12 μL of Negative Control (water)
Your unknown sample reaction	<ul style="list-style-type: none">• Add 18 μL of Premix Solution• Add 10 μL of unknown sample• Add 2 μL of Negative Control (water)
Inhibition-control reaction	<ul style="list-style-type: none">• Add 18 μL of Premix Solution• Add 10 μL of unknown sample• Add 2 μL of the Discriminatory Positive Control (DPC)
Positive control reaction	<ul style="list-style-type: none">• Add 18 μL of Premix Solution• Add 2 μL of the DPC• Add 10 μL of Negative Control (water)

Note: The MycoSEQ™ *Mycoplasma* Discriminatory Positive/Extraction Control can be used as a spike control that is added to the test sample or lysate before sample preparation.

For units:

- With standard 0.2-mL block – Dispense into a standard optical 96-well plate (Cat. No. 4306737).
 - With Fast 0.1-mL block – Dispense into a Fast optical 96-well plate (Cat. No. 4346906).
2. If using a standard 0.2-mL plate, mix each sample by gently pipetting up and down.
 3. Seal the plate with MicroAmp™ Optical Adhesive Film. See “Seal the plates” on page 28.
 4. Briefly centrifuge the reaction plate.



Start the run

1. Double-click Mycoplasma SEQ Example Setup.eds to open the example experiment file you created in “Save the experiment” on page 11.
2. Load the reaction plate into the instrument.
3. To start your 7500 Fast instrument:
 1. Click **Run** in the navigation pane.
 2. Click **START RUN** at the top of any run screen.

Analyze the results


The acceptance criteria that are provided in this section are based on our current knowledge of assay performance in detection of *Mycoplasma* recovered from a wide variety of test sample matrices. We recommend that you qualify and validate the assay internally using samples that are specific to your process and manufacturing environment (raw materials, bioreactor, or cell line samples) to ensure that these criteria are appropriate.

For specific sample types, it may be necessary to make slight changes to the acceptance criteria based on specific results. We can provide you with one-on-one support during this process.

Set the baseline and threshold values

For all reactions, use the default Analysis Settings:

1. Select **Manual C_t**, then set Threshold to **0.2**.
2. Select **Manual Baseline**, then enter the following settings:
 - Start (cycle): **3**
 - End (cycle): **15**





Note: Autobaseline can also be used. To edit the baseline go to **Analysis ▶ Analysis Settings**.

Review the Results Summary

AccuSEQ™ v2.0 software uses the acceptance criteria in Table 3 to provide an automated call summary for each reaction. Use the Call Summary, Plate Layout, and Table views in the **Results Summary** screen to review the experiment results.

1. From the navigation pane, select **Results ▶ Results Summary**.
2. Review the **Call Summary** for results.

Call Summary					
Positive Controls:	0	✓ Pass	0	✗ Fail	
Negative Controls:	0	✓ Pass	0	✗ Fail	
Unknowns:	0	⊕ Present	0	⊖ Absent	0 ? Review 0 ✗ Fail

3. (Optional) Adjust the *Mycoplasma* presence/absence analysis to meet your method qualifications. In the **Call Settings** tab, modify the values of the threshold cycle (C_t), derivative value (DV), and melting temperature (T_m) values. Use the slider interface at the top-half of the screen, to automatically update the comprehensive table at the bottom half of the screen.

Analysis Settings for 18May2017

Call Settings

Unknowns Positive Control and Inhibition Control Negative Control

Review and edit the Mycoplasma Presence Absence call settings for Unknown

C_t threshold for Unknown samples: 30.0 35.0 40.0

DV Range for Unknown samples: 0.2 0.4 0.6 0.8 1.0

T_m Range for Unknown samples: 70.0 75.0 80.0 85.0

Arrow keys can be used to move the slider for finer granularity.

Call	C_t	DV	T_m	Inhibition
Present	< 36.23	≥ 0.8	$75.0 \leq X \leq 82.0$	Pass or Fail or Not Used
Review	< 36.23	$0.8 > X \geq 0.4$	$75.0 \leq X \leq 82.0$	Pass or Fail or Not Used
Absent	< 36.23	< 0.4	$75.0 \leq X \leq 82.0$	Pass
Review	< 36.23	< 0.4	$75.0 \leq X \leq 82.0$	Fail
Absent	< 36.23	< 0.4	$75.0 \leq X \leq 82.0$	Not Used
Absent	≥ 36.23	< 0.4	$75.0 \leq X \leq 82.0$	Pass
Absent	≥ 36.23	< 0.4	$75.0 \leq X \leq 82.0$	Not Used

Revert to Default Admin Analysis Settings
 Revert to Default factory Analysis Settings
 Save as Default Admin Analysis Settings
 Apply Analysis Settings
 Cancel



Guidance for test samples

The table shows criteria for positive and negative calls. A positive call indicates that at least one genome copy of *Mycoplasma* DNA was present in the test reaction and the sample is positive for the presence of *Mycoplasma*. The automated threshold setting for derivative value (DV) of 0.8 for AccuSEQ™ software v2.0 (or later) is equivalent to the 0.05 setting for SDS v1.4 (or later) software.

Note: The values in the tables are subject to your own validation.

Table 2 Recommended acceptance criteria for test samples: AccuSEQ™ software v2.0 or later

Result	C _t	T _m	DV
Positive	< 36.23 C _t	75 – 82°C	≥ 0.8
Negative	≥ 36 C _t	< 75°C	N/A


Guidance for controls



The values in the tables are subject to your own validation.

Table 3 Recommended acceptance criteria for controls: AccuSEQ™ software v2.0 or later.

Control	C _t	T _m	DV
PCR positive control	< 36.23 C _t	82–86°C	> 0.8
Extraction spike control	< 36.23 C _t	82–86°C	> 0.8
No template control	≥ 36.23 C _t	< 75°C	N/A
Blank extraction control	≥ 36.23 C _t	< 75°C	N/A
Inhibition control	ΔC _t < 2	82–86°C	N/A

Guidance for inconclusive results with AccuSEQ™ software v2.0

If a MycoSEQ™ assay does not meet all of the criteria for a positive or negative automatic call, the well displays  (inconclusive). For information about these results:

- Click  (Quality Summary) in the Results navigation pane of the AccuSEQ™ software v2.0 screen.
- Click  (Help) in the toolbar at the top of the AccuSEQ™ software v2.0 screen.
- See Appendix B, “Troubleshooting”.
- Refer to the *AccuSEQ™ Real-Time PCR Detection Software Mycoplasma SEQ Experiments Getting Started Guide*.



Example results with AccuSEQ™ Software v. 2.1.1

Positive control




Figure 1 PCR positive control spiked with 2,000 copies of DPC.



Negative control




Figure 2 No template PCR control.



Blank extraction control




Figure 3 Blank extraction control with PBS.



Positive extraction control




Figure 4 Sample spiked with 2,000 copies of DPC before DNA extraction.



Inhibition control and positive control




Figure 5 No PCR inhibition present; inhibition control and PCR positive control overlaid, with a $\Delta C_t < 2$.



Test sample:
Negative result




Figure 6 Negative result; $C_t > 36.23$ and $T_m < 75^{\circ}\text{C}$.



Test sample:
Positive result




Figure 7 Positive result; $C_t = 15.69$, $T_m > 75^\circ\text{C}$, and Derivative Reporter > 0.8 .



Test sample:
Positive result
with decreased
detection of DPC




Figure 8 Decreased DPC signal can be observed in the presence of very high *mycoplasma* contamination.



PCR inhibition




Figure 9 PCR inhibition, $\Delta C_t > 2$.



Multicomponent plots




Figure 10 An example of a multicomponent plot.




Figure 11 A multicomponent plot with high background signal which will result in PCR inhibition.



Use the kit with 7500 System SDS Software v1.4 or later

Prepare the kit reagents and premix solution

1. Thaw all kit reagents completely.
2. Vortex briefly, then spin down the reagents.
3. Prepare the Premix Solution according to the following table.

Component for premix solution	Volume for one 30- μ L reaction	Volume for four 30- μ L reactions ^[1]
<i>Power</i> SYBR™ Green PCR Master Mix, 2X	15.0 μ L	66.0 μ L
<i>Mycoplasma</i> Real-Time PCR Primer Mix, 10X	3.0 μ L	13.2 μ L
Total premix solution volume	18.0 μL	79.2 μL

^[1] Includes 10% excess to compensate for pipetting errors.

4. Mix the Premix Solution by gently pipetting up and down, then cap the tube.

Prepare the PCR reactions

1. Dispense the following into each well to be used, gently pipetting at the bottom of the well.

To prepare...	In each tube or well...
Negative control reaction	<ul style="list-style-type: none"> • Add 18 μL of Premix Solution • Add 12 μL of Negative Control (water)
Your unknown sample reaction	<ul style="list-style-type: none"> • Add 18 μL of Premix Solution • Add 10 μL of unknown sample • Add 2 μL of Negative Control (water)
Inhibition-control reaction	<ul style="list-style-type: none"> • Add 18 μL of Premix Solution • Add 10 μL of unknown sample • Add 2 μL of the Discriminatory Positive Control (DPC)
Positive control reaction	<ul style="list-style-type: none"> • Add 18 μL of Premix Solution • Add 2 μL of the DPC • Add 10 μL of Negative Control (water)

Note: The MycoSEQ™ *Mycoplasma* Discriminatory Positive/Extraction Control can be used as a spike control that is added to the test sample or lysate before sample preparation

For units:

- With standard 0.2-mL block – Dispense into a standard optical 96-well plate (Cat. No. 4306737).
- With Fast 0.1-mL block – Dispense into a Fast optical 96-well plate (Cat. No. 4346906).

- For each row of wells that you use, place in sequence from left to right the negative control, unknown sample, inhibition control, then positive control. See “Plate layout suggestions” on page 43 for more information.

Pipetting guidelines:

- Use at least one negative and one positive control per run.
- Mix each sample gently by pipetting up and down.
- Use a new tip for each well, even when aliquoting the same solution.

	1	2	3	4	5	6	7	8	9	10	11	12
A	NEG 1 N Mycopl				Sample 1 U Mycopl				Sample 1 IC Mycopl			POS 1 P Mycopl
B					Sample 2 U Mycopl				Sample 2 IC Mycopl			
C					Sample 3 U Mycopl				Sample 3 IC Mycopl			
D					Sample 4 U Mycopl				Sample 4 IC Mycopl			
E					Sample 5 U Mycopl				Sample 5 IC Mycopl			
F					Sample 6 U Mycopl				Sample 6 IC Mycopl			
G					Sample 7 U Mycopl				Sample 7 IC Mycopl			
H					Sample 8 U Mycopl				Sample 8 IC Mycopl			


Wells: U Unknown 8 N Negative Control 1 P Positive Control 1 IC Inhibition Control 8 78 Empty

Figure 12 Example plate layout.


Seal the plates

- Place an optical adhesive cover on the plate, then rub the flat edge of the applicator back and forth along the **long** edge of the plate.

IMPORTANT! Apply significant downward pressure on the applicator to completely seal the wells. Pressure is required to activate the adhesive on the optical cover.



2. Rub the flat edge of the applicator back and forth along the **short** edge (width) of the plate.
3. Rub the edge of the applicator horizontally and vertically between all wells.
4. Rub the edge of the applicator around all outside edges of the plate using small back and forth motions to completely seal around the outside wells.
5. Briefly spin down the plate using a centrifuge with a plate adapter.



IMPORTANT! Make sure that the reagents (and no bubbles) are in the bottom of the wells.

Prepare the plate document

Set up the plate document in the SDS software. For more details, see the *Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Getting Started Guide: Absolute Quantitation using Standard Curve* (Pub. No. 4347825).

1. In the **Assay** drop-down list, select **Absolute Quantification**.
2. Select **SYBR™ detector** with:
 - Quencher Dye set to **none** or **Non Fluorescent**
 - Passive Reference set to **ROX™**
3. Set thermal-cycling conditions as indicated in the table below.

Note: For instruments using the AccuSEQ™ 2.0 Real-Time PCR Software *Mycoplasma* Module, the cycling conditions are pre-programmed in the software.

Step	AmpliTaq Gold™ enzyme activation	PCR		Dissociation ^[1,2,3]			
	HOLD	Cycle (40 cycles)		Melt			
		Denature	Anneal/extend				
Temp	95°C	95°C	60°C	95°C	60°C	95°C	60°C
Time	10 min	15 sec	1 min	15 sec	1 min	15 sec	15 sec

^[1] 7500 and 7500 Fast Systems: from the Instrument tab, click **Add Dissociation Stage** (see Figure 13).

^[2] Applied Biosystems™ Real-Time PCR Instruments: from the Instrument tab, click **Add Dissociation Stage** (see Figure 13). Use default settings..

^[3] For other instruments, refer to their corresponding user guides for dissociation-curve setup information.




Figure 13 The instrument tab for 7500 Fast Real-Time PCR platform with SDS v1.4 software. The run mode is set to Standard 7500.

4. Set Sample Volume to 30 μ L.
5. Select the Standard Run Mode for use with SYBR[™] Green I dye.

Perform PCR

On an Applied Biosystems[™] Real-Time PCR System:

1. Open the plate document that corresponds to the reaction plate ("Prepare the plate document" on page 29).
2. Load the reaction plate into the real-time PCR system.
3. Start the run.

Analyze the results


The acceptance criteria that are provided in this section are based on our current knowledge of assay performance in detection of *Mycoplasma* recovered from a wide variety of test sample matrices. We recommend that you qualify and validate the assay internally using samples that are specific to your process and manufacturing environment (raw materials, bioreactor, or cell line samples) to ensure that these criteria are appropriate.

For specific sample types, it may be necessary to make slight changes to the acceptance criteria based on specific results. We can provide you with one-on-one support during this process.

Set the baseline and threshold values

For all reactions, use the default Analysis Settings:

1. Select **Manual C_t** then set Threshold to **0.2**.
2. Select **Manual Baseline**, then enter the following settings:
 - Start (cycle): **3**
 - End (cycle): **15**



Note: Autobaseline can also be used. To edit the baseline go to **Analysis ► Analysis Settings**.

Guidance for test samples

The table shows criteria for positive and negative calls. A positive call indicates that at least one genome copy of *Mycoplasma* DNA was present in the test reaction and the sample is positive for the presence of *Mycoplasma*. The automated threshold setting for derivative value (DV) of 0.8 for AccuSEQ™ 2.0 (or later) software is equivalent to the 0.05 setting for SDS v1.4 (or later) software.

Note: The values in the tables are subject to your own validation.

Table 4 Recommended acceptance criteria for test samples: SDS software v1.4 or later.

Result	C _t	T _m	DV
Positive	< 36 C _t	75°C – 81°C	≥0.05
Negative	≥ 36 C _t	< 75°C	N/A

Guidance for controls

The values in the tables are subject to your own validation.

Table 5 Recommended acceptance criteria for controls: SDS software v1.4 or later.

Control	C _t	T _m	DV
PCR positive control	< 36 C _t	≈84°C	> 0.05
Extraction spike control	< 36 C _t	≈84°C	> 0.05
No template control	≥ 36 C _t	< 75°C	N/A
Blank extraction control	≥ 36 C _t	< 75°C	N/A
Inhibition control	ΔC _t < 2	≈84°C	N/A

- Both the PCR positive control and the extraction spike control may present extra peaks with T_m < 75°C. These peaks represent primer dimer formation, and they do not interfere with the final results.
- The difference in C_t between the DPC and the inhibition control reaction should be less than 2 cycles. If the unknown sample is negative and the inhibition control shows a ΔC_t > 2 when compared to the positive control, then the PCR is likely inhibited. The sample should be re-purified and the assay repeated. The ΔC_t is calculated by C_t (of inhibition control reaction) – C_t (of positive control reaction).

Example positive results with SDS v1.4 software




Figure 14 *Mycoplasma* contamination (approximately 3×10^6 copies per PCR reaction).




Figure 15 *Mycoplasma* contamination (approximately 2,000 copies per PCR reaction).




Figure 16 *Mycoplasma* contamination (less than 10 copies per PCR reaction).



Example positive control extraction results with SDS v1.4 software




Figure 17 Sample spiked with 2,000 copies of DPC and contaminated with *Mycoplasma* (3×10^6 copies).




Figure 18 Sample spiked with 2,000 copies of DPC and contaminated with *Mycoplasma* (approximately 2,000 copies).




Figure 19 Sample containing 25 copies of *Mycoplasma* and 25 copies of DPC.

Example negative results with SDS v1.4 software




Figure 20 Negative result.





Figure 21 Negative sample spiked with 2,000 copies of DPC.



Troubleshooting

AccuSEQ™ 2.0 software

The following table shows some common reasons for inconclusive results with AccuSEQ™ 2.0 software. For a complete list, click  (Help) in the toolbar at the top of the AccuSEQ™ 2.0 software screen. Refer to the *AccuSEQ™ software: Mycoplasma Getting Started Guide* (Part No. 4425587) for more data analysis information and example results.

Analysis result	Description	Possible cause	Recommended action
Inconclusive	Based on one or more parameter, the software did not make a positive or negative call.	Low sample concentration of <i>Mycoplasma</i> .	<ul style="list-style-type: none">• Call manually according to laboratory guidelines.<i>or</i>• Allow the culture to grow for an additional 24 hours, then repurify the sample and repeat the experiment using assay components that were stored correctly.
Border-line C _t result		Incorrect baseline setting applied	Re-analyze samples using auto baseline, or change to manual start and end baseline settings.
		High SYBR signal from contaminating host cell DNA.	Apply RNase treatment during sample preparation. Contact your local FAS for more information.

Analysis result	Description	Possible cause	Recommended action
Sample inhibits amplification	PCR inhibition shown by: <ul style="list-style-type: none"> Negative unknown sample <i>and</i> Inhibition control $\Delta C_t > 2$ compared to the positive control 	<ul style="list-style-type: none"> Inappropriate sample preparation that results in carryover of chemicals from the lysis buffer. Excess DNA or RNA. SYBR™ signals will have increased signal in Component View. Components in cell culture media or additive (for example, dextran sulfate) may inhibit PCR and change T_m. 	<ul style="list-style-type: none"> Repurify the sample preparation and make sure of appropriate performance of wash and elution that does not carry over chemicals from reagents. Repurify the sample using protocol with RNase and/or DNase treatment. Contact your Field Applications Specialist (FAS) for more information. Most cell culture media inhibit PCR or change T_m without sample preparation. Check if additives to cell culture media inhibit PCR. To address inhibition from Dextran sulfate or Heparin, contact your FAS for more information.
High background signal	High background fluorescence signal: >500,000 fluorescent standard units (FSU).	The sample block is contaminated. Sample may have high concentration of nucleic acid carried from the cell culture during sample preparation.	Run a background calibration to identify the contaminated wells, then decontaminate the sample block. Repeat the experiment using assay components that were stored correctly.

MycoSEQ™ kit

Observation	Possible cause	Action
No positive control or target-specific SYBR™ Green dye signal is detected in inhibition control and/or positive control wells	Improper storage of Power SYBR™ Green PCR Master Mix.	Repeat the assay using properly stored assay components.
	Improper storage of target-specific <i>Mycoplasma</i> Real-Time PCR Primer Mix (10X).	Avoid freezing and thawing assay components. Protect Power SYBR™ Green PCR Master Mix from light.

Observation	Possible cause	Action
No positive control or target-specific SYBR™ Green dye signal is detected in inhibition control and/or positive control wells	Pipetting error (no premix solution added).	Repeat the assay. Make sure to pipet premix solution into all wells.
	Pipetting error (no positive control added).	Repeat the assay. Make sure to pipet positive control into all positive-control wells.
Target-specific signal is detected in negative control wells	Carryover contamination.	<p>Repeat the assay using fresh aliquots of all reagents and clean pipetting equipment.</p> <p>If the negative control continues to show contamination, repeat the assay using a new kit.</p> <p>If the negative control continues to show contamination, contact your Application Specialist.</p>
	High level of nonspecific product formation.	<p>Check the dissociation curve to confirm. Repeat the assay using properly stored assay components.</p> <p>Avoid freezing and thawing assay components. Protect Power SYBR™ Green PCR Master Mix from light.</p>



Background information

Mycoplasmas are the smallest and simplest self-replicating organisms. Their genome sizes range from about 540 kb to 1300 kb, with a G+C content of 23 mol to 41 mol%. Although mycoplasmas are derived from the gram-positive branch of walled eubacteria, their evolution from these walled bacteria resulted in a substantial reduction in genome size and loss of the functions required for synthesis and maintenance of a bacterial cell wall.

Mycoplasmas are a common bacterial contaminant of cell culture samples. Infection is persistent, difficult to detect and diagnose, and very difficult to cure. Mycoplasmas vary in size from 0.2 μm to 0.8 μm , so they can pass through some filters that are used to remove bacteria. Mycoplasma in infected cell cultures can change many cell processes, including altering cell growth rate, inducing morphological changes or cell transformation, and mimicking virus infection. Cell culture in pharmaceutical production must be *Mycoplasma*-free as required by the U.S. Pharmacopoeia and FDA regulatory requirements. Therefore, there is an absolute requirement for routine, periodic testing of possible contamination of all cell cultures used in pharmaceutical manufacturing. Because mycoplasmas grow slowly (the colonies can take up to 3 weeks to develop), traditional culture methods are unacceptable for rapid high-throughput testing. The recently introduced and validated rapid bacterial testing methods that are used in this kit provide for fast *Mycoplasma* screening.



Kit specificity

Sensitivity

The sensitivity of the PCR using this kit is 1 to 10 copies of the target DNA per reaction. Sensitivity of the assay in real culture samples depends on the quality of the sample preparation method that is used. The sample preparation procedure in the *PrepSEQ™ Sample Preparation Kits for Mycoplasma, MMV, and Vesivirus User Guide* (Pub. No. 4465957) allows you to detect:

- 4 to 10 CFU/mL of *Mycoplasma* from 10 mL of cell culture
or
- 4 CFU/mL of *Mycoplasma* from 1 mL of media

Kit specificity

The MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit can detect more than 90 different *Mycoplasma* species, including *Acholeplasma laidlawii* and *Spiroplasma citri*. The kit does not detect other genera or cell-line DNA.

Inclusivity – detectable species

The kit procedure in this protocol is designed to detect over 90 species, including the 14 shown below in the first table. For a complete list of species, contact Technical Support.

Species	Strain/source
<i>Acholeplasma laidlawii</i>	ATCC 23206D
<i>Mycoplasma arginini</i>	ATCC 23838D
<i>Mycoplasma fermentans</i>	ATCC 19989D
<i>Mycoplasma gallisepticum</i>	ATCC 15302
<i>Mycoplasma genitalium</i>	ATCC 33530D
<i>Mycoplasma hominis</i>	ATCC 23114D
<i>Mycoplasma hyorhinis</i>	ATCC 17981D
<i>Mycoplasma hyponeumoniae</i>	ATCC 25095
<i>Mycoplasma orale</i>	ATCC 23714D
<i>Mycoplasma pirum</i>	ATCC 25960D



Species	Strain/source
<i>Mycoplasma pneumoniae</i>	ATCC 15531D
<i>Mycoplasma salivarium</i>	ATCC 23064D
<i>Mycoplasma sinoviae</i>	ATCC 25204
<i>Spiroplasma citri</i>	ATCC 27556D

Exclusivity – undetectable organisms

Organism	Strain/source
<i>Bacillus cereus</i>	ATCC 10876
<i>Bacillus subtilis</i>	ATCC 6051
<i>Campylobacter jejuni</i>	ATCC 29428
<i>Citrobacter freundii</i>	6879
<i>Clostridium perfringens</i>	ATCC 12915
<i>Enterobacter aerogenes</i>	Q87
<i>Enterobacter sakazaki</i>	ATCC 51329
<i>Enterococcus faecalis</i>	ATCC 29212
<i>Escherichia coli</i> O157:H7	43888
<i>Klebsiella oxytoca</i>	ATCC 43165
<i>Lactobacillus bulgaris</i>	ATCC 11842
<i>Listeria ivanovii</i>	ATCC 19119
<i>Listeria monocytogenes</i>	ATCC 7644
<i>Pseudomonas aeruginosa</i>	ATCC 27853
<i>Pseudomonas aeruginosa</i>	ATCC 17423
<i>Shigella</i>	Sfla 395
<i>Shigella</i>	SFL 153
<i>Shigella dysenteriae</i>	ATCC 13313
<i>Shigella dysenteriae</i>	ESCL7-JHH
<i>Staphylococcus aureus</i>	ATCC 43300
<i>Staphylococcus aureus aureus</i>	PE491
<i>Streptococcus faecalis</i>	ATCC 9790



Appendix D Kit specificity
Exclusivity – undetectable organisms

Organism	Strain/source
<i>Vibrio cholerae</i>	036
<i>Yersinia enterocolitica</i>	ATCC 9610
Cat	Novagen™, Cat. No. 69235-3
Cow	Novagen™, Cat. No. 69238-3
Chicken	Novagen™, Cat. No. 69233-3
Chimpanzee	Bios, Inc. ^[1]
CHO	ATCC CCL-61
HeLa	ATCC CCL-2
Horse	Pel-Freez Biologicals, Cat. No. 39339-5
Orangutang	Bios, Inc. ^[1]
Pig	Novagen™, Cat. No. 69230-3
Rabbit	Pel-Freez Biologicals, Cat. No. 31130-1
Rat	Novagen™, Cat. No. 69238-3
Sheep	Novagen™, Cat. No. 69231-3

^[1] No longer available



Good PCR practices

PCR assays require special laboratory practices to avoid false positive amplifications. The high throughput and repetition of these assays can lead to amplification of one DNA molecule. Follow the guidelines below to prevent contamination and nonspecific amplification.

Good laboratory practices for PCR and RT-PCR

When preparing samples for PCR or RT-PCR amplification:

- Wear clean gloves and a clean lab coat.
 - Do not wear the same gloves and lab coat that you have previously used when handling amplified products or preparing samples.
- Change gloves if you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation and reaction setup.
 - Amplification and analysis of products.
- Do not bring amplified products into the reaction setup area.
- Open and close all sample tubes carefully. Avoid splashing or spraying samples.
- Keep reactions and components capped as much as possible.
- Use a positive-displacement pipettor or aerosol-resistant barrier pipette tips.
- Clean lab benches and equipment periodically with 10% bleach solution or DNA decontamination solution.

Plate layout suggestions

- For each plate row, dispense in sequence from left to right: negative controls, unknown samples, inhibition controls, and positive controls (at the end of the row or column).
- Place positive controls in one of the outer columns.
- If possible, separate all samples from each other by at least one well. If space is limited, place at least one well between unknown samples and controls.

Documentation and support

Related documentation

Document	Pub. No.	Description
<i>MycoSEQ™ Mycoplasma Real-Time PCR Detection Kit Quick Reference</i>	4465876	Provides brief, concise instructions on using the MycoSEQ™ Mycoplasma Detection Kit.
<i>ViralSEQ™ Mouse Minute Virus (MMV) Real-Time PCR Detection Kit Quick Reference</i>	4445236	Provides brief, concise instructions on using the ViralSEQ™ Mouse Minute Virus Real-Time PCR Detection Kit.
<i>ViralSEQ™ Mouse Minute Virus (MMV) Real-Time PCR Detection Kit User Guide</i>	4445235	Describes the ViralSEQ™ Mouse Minute Virus Real-Time PCR Detection Kit and provides information on preparing, running, and troubleshooting MMV detection.
<i>PrepSEQ™ Sample Preparation Kits for Mycoplasma, MMV, and Vesivirus Quick Reference</i>	4465875	Provides brief, concise instructions on using the PrepSEQ™ Sample Preparation Kits.
<i>PrepSEQ™ Sample Preparation Kits for Mycoplasma, MMV, and Vesivirus User Guide</i>	4465957	Describes the PrepSEQ™ Sample Preparation Kits and provides information on preparing, running, and troubleshooting sample preparation.
<i>PrepSEQ™ Nucleic Acid Extraction Kit Quick Reference</i>	4406303	Provides brief, concise instructions on using the PrepSEQ™ Nucleic Acid Extraction Kit.
<i>PrepSEQ™ Nucleic Acid Extraction Kit User Guide</i>	4400739	Describes the PrepSEQ™ Nucleic Acid Extraction Kit and provides information on preparing, running, and troubleshooting nucleic acid extractions.
<i>Introduction to TaqMan® and SYBR™ Green Chemistries for Real-Time PCR Protocol</i>	4407003	Describes the TaqMan® and SYBR™ Green Chemistries for Real-Time PCR and provides information on preparing, running, and troubleshooting PCR.

Document	Pub. No.	Description
<i>AccuSEQ™ software: Mycoplasma Getting Started Guide</i>	4425587	Provides brief, step-by-step procedures for <i>Mycoplasma</i> detection. It is designed to help you quickly learn to use the AccuSEQ™ Real-Time PCR Detection Software for Mycoplasma SEQ Experiments.
<i>Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Getting Started Guide: Absolute Quantitation using Standard Curve</i>	4347825	Provides brief, step-by-step procedures for absolute quantitation using a standard curve. It is designed to help you quickly learn to use the Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System.

For information on new assays and updated product documentation, go to **thermofisher.com**.

Portable document format (PDF) versions of this guide and the documents listed above are available at **thermofisher.com**

Customer and technical support

Visit **thermofisher.com/support** for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.thermofisher.com/us/en/home/global/terms-and-conditions.html**. If you have any questions, please contact Life Technologies at **www.thermofisher.com/support**.

