

ProteinSEQ™ Protein A Quantification Kit

USER GUIDE

Catalog Number 4469343

Publication Number MAN0013524

Revision B.0

The information in this guide is subject to change without notice.

DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity

Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

Trademarks

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Adobe and Reader are registered trademarks of Adobe Systems Incorporated. Microsoft and Excel are registered trademarks of Microsoft Corporation. GraphPad Prism is a registered trademark of GraphPad Software, Inc. MabSelect SuRe is a trademark of GE Healthcare.

©2015 Thermo Fisher Scientific Inc. All rights reserved.

Contents

About this guide	5
Revision history	5
■ Product information	6
Product description	6
Kit contents and storage	6
Materials and equipment not provided	7
■ Methods	9
Workflow	9
Important procedural guidelines	9
Important serial dilution guidelines	10
Before you begin (first time use only)	10
Before you begin each time	11
Pre-dilute the Protein A standard	11
Pre-dilute the Repligen Protein A standard	11
Pre-dilute the GE Healthcare MabSelect SuRe™ Protein A standard	11
Prepare serial dilutions of the diluted Protein A standard	12
Prepare samples	13
Adjust sample pH and/or concentration	13
Treat samples with 2X Protein A Sample Prep Reagent	13
Prepare plates for the MagMAX™ Express-96 run	14
Label plates	14
Prepare wash plates	15
Prepare qPCR plate	15
Prepare probes plate	15
Prepare the capture plate	16
Run plates in the MagMAX™ Express-96 Magnetic Particle Processor	17
Run qPCR reaction	19
Set up and run qPCR on the 7500 Fast instrument with AccuSEQ™ software v2.0 or later	20
Perform data analysis	26
Perform data analysis with AccuSEQ™ software v2.0 or later	26
Perform data analysis without AccuSEQ™ software	26

■	APPENDIX A	Troubleshooting	28
■	APPENDIX B	Design guidelines for ProteinSEQ™ System	
		Protein A spike experiments	30
		About spike experiments	30
		Important experimental design considerations	32
		Guidelines for spike input concentration	32
		Guidelines for matrix selection	33
		Protein A quantitation example	34
■	APPENDIX C	Good laboratory practices for PCR and RT-PCR	36
■	APPENDIX D	Safety	37
		Chemical safety	38
		Biological hazard safety	39
		Documentation and support	40
		Related documentation	40
		Customer and technical support	40
		Limited product warranty	41

About this guide

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Revision history

Revision	Date	Description
A.0	September 2015	New user guide.
B.0	December 2015	Corrected qPCR well volume setting for standard PCR plates. Added “v2.0 or later” to AccuSEQ™ software references



Product information

Product description

Use the ProteinSEQ™ Protein A Quantification Kit to quantify the amount of leached Protein A present in your bioprocess sample(s). The ProteinSEQ™ Protein A workflow consists of:

- Plate preparation
- A semi-automated sample processing run on the MagMAX™ Express-96 Magnetic Particle Processor
- A qPCR run on an Applied Biosystems™ 7500 Fast Real-Time PCR System (or equivalent system)
- Data analysis using AccuSEQ™ Real-Time PCR software v2.0 or later, or equivalent curve fitting software

Kit contents and storage

Table 1 ProteinSEQ™ Protein A Quantification Kit (Cat. no. 4469343)

Box Part No.	Contents	Cap color	Item Part No.	Quantity	Volume	Storage conditions
A25001	Repligen Protein A Standard (1 mg/mL)	●	4469347	1	0.1 mL	-20°C
	Protein A 5' Probe	●	4469376	1	0.7 mL	
	Protein A 3' Probe	●	4469381	1	0.7 mL	
	2X Protein A Sample Prep Reagent	○	4485137	1	10 mL	
	ProteinSEQ™ Ligation and Assay Mix	●	100024250	1	0.6 mL	
	ProteinSEQ™ Ligase	●	A25020	1	55 µL	
A25002	Fast Master Mix, 2X	○	4448622	1	6 mL	4°C
	Wash Buffer	○	4469398	1	100 mL	
	Protein A Capture Beads	○	4469368	1	1.2 mL	
	ProteinSEQ™ Elution Buffer	○	100024368	1	7 mL	
	ProteinSEQ™ Diluent	○	4469405	1	41 mL	



Materials and equipment not provided

Unless otherwise indicated, all materials are available through **thermofisher.com**.

MLS: Fisher Scientific (www.fisherscientific.com) or other major laboratory supplier.

Item	Source
Instruments	
Applied Biosystems™ 7500 Fast Real-Time PCR System	4365464
Equipment	
MagMAX™ Express-96 Deep Well Magnetic Particle Processor	4456933
MagMAX™ Express-96 PCR Well Magnetic Head	4472991
qPCR Plate Adaptor for MagMAX™ Express-96 Processor	Contact Technical Support
Heat block for 1.5/2-mL tubes	MLS
Benchtop microcentrifuge	MLS
Plate centrifuge	MLS
Benchtop vortexer	MLS
Reagents	
MabSelect SuRe™ Protein A, 2 mg/mL ^[1]	GE Healthcare 28-4018-60 (custom order)
0.5 M Sodium Phosphate Buffer, pH 7	MLS
Consumables	
Aerosol-resistant pipette tips	MLS
MagMAX™ Express Tip Comb	4472784
Disposable gloves	MLS
25-mL reagent reservoir	VistaLab Technologies 3054-1002 or equivalent
15-mL conical tube	AM12500 or equivalent
Pipettors <ul style="list-style-type: none"> • P20, P200, and P1000 single-channel pipettors • P200 and P1000 multichannel pipettors, 8- or 12-channel 	MLS
qPCR plates <ul style="list-style-type: none"> • MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode [0.1 mL]^[2] or • MicroAmp™ Optical 96-Well Reaction Plate with Barcode [0.2 mL]^[3] 	<ul style="list-style-type: none"> • 4346906 or • 4306737




Item	Source
Consumables (continued)	
Additional plates <ul style="list-style-type: none">MagMAX™ Express-96 Skirted Low Profile PlatesMagMAX™ Express-96 Standard Plates	<ul style="list-style-type: none">4472783 (also available as AB-0800 from www.fisherscientific.com)4388475
MicroAmp™ Optical Adhesive Film	4360954
Nonstick, RNase-free Microfuge Tubes, 1.5 mL	AM12450
Sterile Microcentrifuge Tubes with Screw Caps	Fisher Scientific 02-681-371 or equivalent
Aluminum foil for heat block	MLS
Software	
AccuSEQ™ Real-Time PCR software v2.0 or later, or equivalent curve fitting software	4443420 or MLS
Microsoft™ Excel™ software	www.microsoft.com

- [1] Optional, if not using Protein A provided in kit.
[2] For use with thermal cycler FAST sample blocks.
[3] For use with thermal cycler standard sample blocks.



Methods

Workflow



Important procedural guidelines

- The MagMAX™ Express-96 PCR Well Magnetic Head is fragile. Handle with care.
- Use serially diluted standards when performing spiking studies. See Appendix B, “Design guidelines for ProteinSEQ™ System Protein A spike experiments” for spiking guidelines.
- Run all reactions in triplicate.
- We recommend digital multi-channel pipettors for transfers into the MagMAX™ Express plates.
- Working solutions and plates can be kept at room temperature during assay setup.



Important serial dilution guidelines

Pre-dilution guidelines— Pre-dilute the Protein A standard to obtain a concentration (40 ng/mL) to use as the starting point for the dilution series.

- Prepare the pre-dilutions in 1.5-mL non-stick RNase-free microfuge tubes (Cat. no. AM12450 or equivalent).
- If you are using:
 - Repligen Protein A Standard (1 mg/mL) provided with the ProteinSEQ™ Protein A Quantification Kit, see “Pre-dilute the Repligen Protein A standard” on page 11.
 - GE Healthcare MabSelect SuRe™ Protein A (2 mg/mL; GE Healthcare 28-4018-60, see “Pre-dilute the GE Healthcare MabSelect SuRe™ Protein A standard” on page 11.

Serial dilution guidelines — Prepare serial dilutions using the pre-diluted Protein A standard (D4).

- Prepare serial dilutions in a MagMAX™ Express-96 Standard Plate (Cat. no. 4388475 or equivalent polypropylene 96-well plate). If your standard dilution volume is >300 µL, prepare the serial dilutions in 1.5-mL non-stick RNase-free microfuge tubes (Cat. no. AM12450 or equivalent).
- Use a new pipette tip for each transfer.
- Pipet gently to minimize foaming and/or bubble formation.
- It is critical to mix standards during serial dilution. After each transfer,
 - **If preparing serial dilutions in microfuge tubes** — Invert the tube several times to mix.
 - **If preparing serial dilutions in a 96-well plate** — Gently pipet up and down 5-8 times to increase mixing efficiency.

Before you begin (first time use only)

- Contact your local MagMAX™ sales or service representative to prepare your MagMAX™ Express-96 Magnetic Particle Processor for use with ProteinSEQ™ assays and to obtain the following items:
 - The ProteinSEQ™ Protein A script for the MagMAX™ Express-96 Magnetic Particle Processor (upload before you perform a ProteinSEQ™ assay for the first time)
 - The appropriate PCR Plate Adaptor (Fast or Standard)
 - If you are using GraphPad™ software for data analysis, the Protein A Master Template
- Ask your local representative if your MagMAX™ Express-96 Magnetic Particle Processor supports plate hold-downs. If supported, your local representative should install the plate hold-downs before you perform a ProteinSEQ™ assay for the first time.
- Calibrate your heat source. It is critical to heat samples to 100°C during incubation. Experiments performed at high altitude may require further assay development such as extended heating time or multiple heating and cooling cycles.



Before you begin each time

- Thaw the 2X Protein A Sample Prep Reagent at room temperature for an hour before preparing samples.
- Before preparing the plates for each assay, clean the pipettors, plate racks and the microcentrifuge (if using tubes for standard dilution) to avoid cross-contamination.

Pre-dilute the Protein A standard

Use the appropriate procedure for your Protein A standard.

Note: See “Important serial dilution guidelines” on page 10.

Pre-dilute the Repligen Protein A standard

If you are using the Repligen Protein A Standard provided with the kit, pre-dilute the standard from 1 mg/mL to 40 ng/mL:

1. Dispense the following to four labeled 1.5-mL microfuge tubes:

Dilution Tubes ^[1]	ProteinSEQ™ diluent	Repligen Protein A standard	Final volume	Final Protein A concentration
D1	480 µL	20 µL (1 mg/mL)	500 µL	40 µg/mL
D2	900 µL	100 µL from D1	1000 µL	4 µg/mL
D3	900 µL	100 µL from D2	1000 µL	0.4 µg/mL
D4	900 µL	100 µL from D3	1000 µL	40 ng/mL

^[1] The pre-dilutions (D1 through D4) can be frozen at -20°C and thawed up to 5 times over one month.

2. After you add the Protein A standard to each tube, vortex the tube for 2 seconds, then briefly spin before making the next transfer.

Pre-dilute the GE Healthcare MabSelect SuRe™ Protein A standard

If you are using GE Healthcare MabSelect SuRe™ Protein A standard (not provided with kit), pre-dilute the standard from 2 mg/mL to 40 ng/mL:

1. Dispense the following to four labeled 1.5-mL microfuge tubes:

Dilution Tubes	ProteinSEQ™ diluent	MabSelect SuRe™ Protein A Standard	Final volume	Final Protein A concentration
D1	490 µL	10 µL (2 mg/mL)	500 µL	40 µg/mL
D2	900 µL	100 µL from D1	1000 µL	4 µg/mL
D3	900 µL	100 µL from D2	1000 µL	0.4 µg/mL
D4	900 µL	100 µL from D3	1000 µL	40 ng/mL

2. After you add the Protein A standard to each tube, vortex the tube for 2 seconds, then briefly spin before making the next transfer.



Prepare serial dilutions of the diluted Protein A standard

See “Important serial dilution guidelines” on page 10.

Prepare a dilution series using D4 from the previous step.

1. Label a MagMAX™ Express-96 Standard Plate “Dilution”.
2. Aliquot approximately 2 mL of ProteinSEQ™ Diluent into a 25-mL reagent reservoir.
3. Dispense 200 μ L ProteinSEQ™ Diluent to wells A1 through H1 of the labeled plate. Use a multi-channel pipette.
4. Add 50 μ L diluted Protein A Standard (40 ng/mL; D4 from “Pre-dilute the Protein A standard” on page 11) to the A1 plate well, then pipet up and down 5 times to mix.
5. Transfer the remaining dilutions. Use a new pipette tip for each transfer. After each transfer, pipet up and down 5 times to mix. See Table 2.





Table 2 Prepare serial dilutions

Serial dilutions	ProteinSEQ™ Diluent	Dilution transfer	Protein A concentration
SD1	200 μ L	50 μ L of diluted Protein A Standard (40 ng/mL)	8000 pg/mL
SD2	200 μ L	50 μ L from SD1	1600 pg/mL
SD3	200 μ L	50 μ L from SD2	320 pg/mL
SD4	200 μ L	50 μ L from SD3	64 pg/mL
SD5	200 μ L	50 μ L from SD4	12.8 pg/mL
SD6	200 μ L	50 μ L from SD5	2.56 pg/mL
SD7	200 μ L	50 μ L from SD6	0.51 pg/mL
NPC	200 μ L	0	0



Prepare samples

The following procedures apply to the unknown samples and do not need to be performed for the Protein A standards used to generate the standard curve.

Adjust sample pH and/or concentration

For best results, adjust the sample pH, IgG concentration, and buffer concentration according to the table:

Sample type	Recommendation
pH < 6	Adjust the pH to > 6 with 0.5M Sodium Phosphate buffer pH 7
pH > 8	Adjust the pH to <8 with 0.5M Sodium Phosphate buffer pH 7
IgG > 25 mg/mL	Dilute to < 25 mg/mL with ProteinSEQ™ Diluent
Carboxylate buffer (Citrate)	Dilute to < 50 mM with ProteinSEQ™ Diluent

IMPORTANT! Do not dilute samples with 2X Protein A Sample Prep Reagent; the final concentration of the 2X Protein A Sample Prep Reagent must be 1X. If necessary to dilute samples, use the ProteinSEQ™ Diluent.

Treat samples with 2X Protein A Sample Prep Reagent

Heating the samples causes the IgG to denature and precipitate. The IgG pellets during centrifugation. The Protein A molecules remain in the supernatant.

Note: Use screw-cap tubes during heat incubation.

1. In a 1.5-mL screw cap or safe-lock tube, add 150 µL of sample to 150 µL of 2X Protein A Sample Prep Reagent.

Note: 2X Protein A Sample Prep Reagent and ProteinSEQ™ Elution Buffer bottle are the same size and cap color. Confirm that you are using the correct reagent.

2. Vortex for 5–10 seconds, then quick spin to collect the contents at the bottom of the tube.
3. Place the sample tubes in a heat block, cover the tubes with aluminum foil, then incubate the tubes at 100°C for 10 minutes.




CAUTION! Adding water to heat blocks may create an electrical hazard. Refer to the manufacturer's guidelines for use of water in heat blocks.

4. Remove the tubes from the heat block and cool at room temperature for 5 minutes.
5. Vortex the tubes on a bench-top vortexer for 3 seconds at maximum speed.
6. Centrifuge the tubes for 5 minutes at 16,000 × g (~14,000 rpm).



- Transfer the entire volume of supernatant to an empty well in the Dilution plate that contains the serially diluted standards.







Prepare plates for the MagMAX™ Express-96 run

Note: When preparing wash, capture, qPCR, and probes plates, dispense at bottom of wells to prevent bubble formation (bubbles prevent effective mixing during the MagMAX™ Express-96 run). If bubbles form, quick-spin the plate at $560 \times g$ (~2000 rpm) in a plate centrifuge.

Label plates

- Label 8 plates:


No. of plates	Plate type		Cat. no.	Label(s)
6	MagMAX™ Express-96 Skirted Low Profile Plates		4472783	Capture Probes Wash 1 Wash 2 Wash 3 Wash 4
1	MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode [0.1 mL] ^[1] or MicroAmp™ Optical 96-Well Reaction Plate with Barcode [0.2 mL] ^[2]	 	4346906 or 4306737	qPCR
1	MagMAX™ Express-96 Standard Plate (200 µL)		4388475 ^[3]	Comb

^[1] For use with thermal cycler FAST sample blocks; shown in Fast PCR Plate Adaptor.

^[2] For use with thermal cycler standard sample blocks; shown in Standard PCR Plate Adaptor.


^[3] Or equivalent polypropylene plate.

- Insert the plate labelled “qPCR” into the appropriate PCR Plate Adaptor [Fast (on left) or Standard (on right); request from your local sales or service representative].






- Place a MagMAX™ Express PCR Head Tip Comb (Cat. no. 4472784) in the plate labelled “Comb”.



Prepare wash plates

- Pour approximately 30 mL Wash Buffer into a fresh reagent reservoir.
- Dispense 100 µL Wash Buffer into each well of the 4 wash plates with a multi-channel pipette.



Prepare qPCR plate




Dispense ProteinSEQ™ Elution Buffer into each well of the qPCR plate:

- Fast PCR plate — 15 µL per well
- Standard PCR plate — 25 µL per well

Note: 2X Protein A Sample Prep Reagent and ProteinSEQ™ Elution Buffer bottle are the same size and cap color. Confirm that you are using the correct reagent.


Prepare probes plate

- Add the assay probe reagents to a 15-mL tube in the order shown in the table. Scale the volumes as needed for the number of reactions, including recommended overages. Vortex for 3 seconds at medium speed, then keep the 15-mL tube on ice.

Reagent	Cap color	Volume ^[1]		
		1 rxn	48 rxn	96 rxn
ProteinSEQ™ Diluent	Clear 	56.3 µL	2700 µL	5400 µL
Protein A 5' Probe	Grey 	3.1 µL	150 µL	300 µL
Protein A 3' Probe	Yellow 	3.1 µL	150 µL	300 µL
Total		62.5 µL	3000 µL	6000 µL

^[1] Includes 25% overage.

- Invert the assay probe mix tube several times to mix, transfer to a reagent reservoir, then dispense 50 µL assay probe mix into each well of the Probes plate with a multi-channel pipette.





Prepare the capture plate

1. Vortex the Protein A Capture Beads for 3 seconds at medium speed 2 times.
2. In a 15-mL conical tube, mix the Protein A Capture Beads and the ProteinSEQ™ Diluent using the volumes from the table:

Reagent	Cap color	Volume ^[1]		
		1 rxn	48 rxn	96 rxn
ProteinSEQ™ Diluent	Clear	21.3 µL	1020 µL	2040 µL
Protein A Capture Beads	Clear	3.8 µL	180 µL	360 µL
Total		25.1 µL	1200 µL	2400 µL

^[1] Includes 25% overage.

3. Immediately dispense 20 µL of diluted Protein A Capture Beads into each well of the Capture plate. Use a multi-channel pipette.
4. Transfer 50 µL of each standard and sample from the Dilution plate to the capture plate in triplicate. Use a multi-channel pipette.
The final volume in the capture plate is 70 µL per well.





Figure 1 Transfer standards and samples to the capture plate

Example: Add 50 µL of each Protein A standard dilution and sample from the Dilution plate to the Protein A Capture Beads in the Capture plate to create 3 replicates of each standard dilution and sample.



Run plates in the MagMAX™ Express-96 Magnetic Particle Processor

1. Turn on the MagMAX™ Express-96 instrument and select the **Protein A** program from the screen.

Note: The instrument automatically resets each time you turn it on.

2. Press **START** to initiate plate loading. Follow the prompts on the display screen to load each plate onto the MagMAX™ Express-96 turntable, starting with “Comb” (see Figure 2). Slide each plate into the plate hold-down (if present).

IMPORTANT! When loading the **Tip Comb** in position 8, confirm that it rests in a MagMAX™ Express 96-well Standard Plate (200 µL; Cat. no. 4388475), *not* a MagMAX™ Express-96 Skirted Low Profile Plate.

For all plates, verify that A1 on the plate aligns with A1 on the instrument.




Figure 2 Plate positions in the MagMAX™ Express-96 turntable




Figure 3 Loading plate (with plate hold-down)



3. Load the last (Capture) plate, then press **START** to begin the run.
The run requires ~2 hours. When the run is complete, the screen displays “Proceed to qPCR”.

Plate	Step	Time
Capture	Protein A binds to Capture Beads	65 min
Wash 1	Capture Beads are washed	4.5 min
Wash 2	Capture Beads are washed	4.5 min
Probes	Probe binds to Protein A on Capture Beads	35 min
Wash 3	Capture Beads are washed	4.5 min
Wash 4	Capture Beads are washed	4.5 min
qPCR	Beads are released into qPCR plate	20 sec
Total		120 min

4. When the MagMAX™ Express-96 program is complete, carefully remove the qPCR plate. Discard the Capture, Wash, and Probes plates.

IMPORTANT! Do not discard the plate adaptor.

Note: Discard the PCR Head Tip Comb. The plate used to hold the PCR Head Tip Comb can be reused in future runs.

The qPCR plate contains Protein A Capture Beads in ProteinSEQ™ Elution Buffer (total volume 30 µL).

Proceed immediately to “Run qPCR reaction” on page 19.






Run qPCR reaction

IMPORTANT! ProteinSEQ™ detection is based on qPCR, which is a highly sensitive technique with potential for cross-contamination. After the run completes, dispose of the qPCR plate. Do not remove the optical film from the qPCR plate; removing the film introduces amplicon contamination into the local environment. See Appendix C, “Good laboratory practices for PCR and RT-PCR”.

1. Prepare the Ligation/qPCR mix in a 15-mL tube according to the volumes shown in the appropriate table, including recommended overages, then briefly vortex to mix.




Table 3 Reagent volumes for **FAST** PCR plates

Reagent	Cap color	Volumes ^[1]		
		1 rxn	48 rxn	96 rxn
Fast Master Mix, 2X	Clear 	20 µL	960 µL	1920 µL
ProteinSEQ™ Ligation and Assay Mix ^[2]	Green 	2 µL	96 µL	192 µL
ProteinSEQ™ Ligase	Orange 	0.2 µL	9.6 µL	19.2 µL
Total		22.2 µL	1065.6 µL	2131.2 µL

^[1] Includes 35% overage.

^[2] Contains FAM™ dye and primers.

Table 4 Reagent volumes for Standard (**non-FAST**) PCR plates

Reagent	Cap color	Volumes ^[1]		
		1 rxn	48 rxn	96 rxn
Fast Master Mix, 2X	Clear 	32.5 µL	1560 µL	3120 µL
ProteinSEQ™ Ligation and Assay Mix ^[2]	Green 	3.25 µL	156 µL	312 µL
ProteinSEQ™ Ligase	Orange 	0.26 µL	12.5 µL	25 µL
Total		36.01 µL	1728.5 µL	3457 µL

^[1] Includes 30% overage; volumes for 48 and 96 reactions are rounded to nearest tenth.

^[2] Contains FAM™ dye and primers.

2. Transfer the Ligation/qPCR mix to each bead-containing well of the qPCR plate that was prepared on the MagMAX™ Express-96 instrument.
 - **For Fast PCR plates** — Use 15 µL per well
 - **For Standard PCR plates** — Use 25 µL per well

Note: Dispense the mix to the sides of the well; do not mix after dispensing.



3. Seal the qPCR plate with an optical film, centrifuge for 3 seconds at 500 rpm, then load the plate on a 7500 Fast Real-Time PCR System (or equivalent).
4. Set up the qPCR run.
 - If you are using the AccuSEQ™ system software, follow the setup instruction in “Set up and run qPCR on the 7500 Fast instrument with AccuSEQ™ software v2.0 or later” on page 20.
 - If you are using other equivalent software (for example, SDS 1.4 software) use the following settings:

–	Stage	Temperature	Time
	Hold	37°C	10 minutes
	Hold	95°C	20 seconds
	40 cycles	95°C	3 seconds
		60°C	30 seconds
–	Setting	FAST plates	Standard (non-FAST) plates
	Protein A standards and sample wells volume	30 µL	50 µL
	Detection dye set	FAM	FAM
	Quencher	none	none


5. Start the run.
6. After the run completes, dispose of the qPCR plate. Do not remove the optical film from the qPCR plate; removing the film introduces amplicon contamination into the local environment.

Set up and run qPCR on the 7500 Fast instrument with AccuSEQ™ software v2.0 or later

1. From the home screen click **Create Custom Experiment**.
2. Make the following selections in the Experiment Properties Pane:
 - a. 7500 Fast (96 Wells)
 - b. Quantitation – Standard Curve
 - c. TaqMan™ Reagents




- d. Fast (~ 50 minutes to complete a run; 10 minutes for ligation and 40 minutes for qPCR)




3. Define Sample Number and Name.

- Click **Plate Setup** in the Experiment Menu Pane.
- Enter **Protein A** as the target name, select **FAM** as the reporter and **None** as the quencher.
- Enter the number and name of your samples, excluding replicates. Click **Add New Sample** to enter the number of samples to be run.
For example, if you have four samples run in triplicate, you would define four samples in this step. Replicates of those four samples will be defined in the next step.





4. Access the Standard Curve Dialog as follows:
 - a. Select the **Assign Targets and Samples** tab.
 - b. Click **Define and Setup Standards**.



5. Setup the Standard Curve as follows:
 - a. In the Define the standard curve tab, enter 7 for “# of Points”, 3 for “# of Replicates”, 8000 for “Starting Quantity”, and 1:5 for “Serial Factor”.
 - b. Click **Let Me Select Wells**. Click, hold, and drag the plate map to select the wells to be used as standards.
 - c. Click **Apply**.




d. Click **Close**.

6. Assign sample name to wells.

- a. In Plate Layout, select all wells that will be assigned as replicates for Sample 1. In this example, A4, A5 and A6 are selected.




- b. Click **Assign** next to the appropriate sample. Repeat for all unknown samples.




In this example, “Sample 1” is assigned to wells A4, A5 and A6.

7. Assign Unknown Well Type as follows:
 - a. Select all wells that will be designated as unknowns.
 - b. Click **Assign** under Assign Targets to the Selected Wells.
 - c. Click the blue **U** to assign the wells as unknowns.
 - d. Verify that **ROX** is selected as the passive reference.






8. Setup qPCR parameters as follows:
 - a. Click **Run Method**.
 - b. Enter a reaction volume of **30 μ L**.
 - c. Verify reaction parameters match those shown:



9. Save the setup as follows:
 - a. Select **File** ► **Save As...** (or select **File** ► **Save As Template...** to save this experimental setup for future use).
 - b. Click **Start Run**.





Perform data analysis

Perform data analysis with AccuSEQ™ software v2.0 or later

1. In the AccuSEQ™ software, select autobaseline **on** and set the threshold manually to **0.2**.
2. Use the AccuSEQ™ software to fit standards to a curve using a non-linear method and obtain interpolated values for the unknowns.
 - 4PL is commonly used for symmetric curves with asymptotes for both the lower and upper Protein A concentrations.
 - 5PL is commonly used if the curve is asymmetric or if either the lower or upper asymptote is not present.
 - Apply 1/Y or 1/Y² weighting according to your criteria.
3. Export the data to a Microsoft™ Excel™ spreadsheet for custom statistical analysis.
4. Evaluate the dynamic range using %CV and the quality of the curve fit.
 - R² is appropriate for judging linear fits but it is not an appropriate metric for evaluating the quality of a non-linear fit.
 - Common acceptance criteria for non-linear curve fits are back-calculation values of 80–120% throughout the curve and 75–125% at the LLOQ.
 - Common acceptance criteria for precision are %CV ≤20% throughout the curve and ≤25% at the LLOQ.
5. Obtain the final concentration for each sample by correcting for sample dilution and spike concentrations, if used.

Perform data analysis without AccuSEQ™ software

1. Select autobaseline **on** and set the C_t threshold manually to **0.2**. Determine the C_t values.
2. Export the raw data from the qPCR software to a Microsoft™ Excel™ spreadsheet, then export from Microsoft™ Excel™ to your fitting program of choice. Transform the values to logarithmic values.

Note: If you use GraphPad™, the Protein A Master Template (a Microsoft™ Excel™ template available from your local sales or service representative) facilitates this process.

3. Fit standards to a curve using a non-linear method and obtain interpolated values for the unknowns.
 - 4PL is commonly used for symmetric curves with asymptotes for both the lower and upper Protein A concentrations.
 - 5PL is commonly used if the curve is asymmetric or if either the lower or upper asymptote is not present.
 - Apply 1/Y or 1/Y² weighting according to your criteria.
4. Transform concentration values from logarithmic to linear values.



5. Evaluate the dynamic range using %CV and the quality of the curve fit.
 - R^2 is appropriate for judging linear fits but it is not an appropriate metric for evaluating the quality of a non-linear fit.
 - Common acceptance criteria for non-linear curve fits are back-calculation values of 80–120% throughout the curve and 75–125% at the LLOQ.
 - Common acceptance criteria for precision are %CV $\leq 20\%$ throughout the curve and $\leq 25\%$ at the LLOQ.
6. Obtain the final concentration for each sample by correcting for sample dilution and spike concentrations, if used.



Troubleshooting

Observation	Possible cause	Recommended action
Capture beads remain on comb	Misalignment of MagMAX™ magnetic head.	Contact your local Technical Support for re-alignment of the instrument.
	qPCR plate not placed in appropriate PCR Plate Adaptor during MagMAX™ Express-96 instrument run.	Use the appropriate PCR Plate Adaptor (Fast or Standard; request from your local sales or service representative).
The standard curve plateaus at the lower standard concentrations and the NPC C_T is less than 28	Cross-contamination of Protein A or ligated product.	Decontaminate the bench and pipettors. Change gloves frequently and follow other good PCR practices. After the run completes, dispose of the qPCR plate. Do not remove the optical film from the qPCR plate; removing the film introduces amplicon contamination into the local environment. See Appendix C, "Good laboratory practices for PCR and RT-PCR". Before preparing the plates for each assay, clean the pipettors, plate racks and the microcentrifuge (if using tubes for standard dilution) to avoid cross-contamination.
		If prone to contamination, change the order of standards, sample, and plate preparation as follows: <ul style="list-style-type: none"> • Label plates • Prepare wash plates • Prepare probes plate • Move the prepared plates near the MagMAX™ instrument. • Prepare standards and samples. • Prepare capture plate
	The reagents are contaminated.	Use new reagents.
The C_T at 1600 pg/mL is above 20 and the NPC C_T is undetermined	Expired kit.	Check kit expiration date.
	Errors in reaction or run setup.	Repeat assay preparation. Make sure that the components are added in the recommended order.
Trending increase in C_T value for standard concentrations from run to run	Deterioration of standards.	Prepare fresh standards. Verify kit expiration date.
Random decrease in C_T during run	Cross-contamination of concentrated standards or samples with lower concentration samples.	Repeat experiment.



Observation	Possible cause	Recommended action
Random failures across the plate	Air bubbles introduced into plate wells during plate setup.	Dispense at bottom of wells to prevent bubble formation (bubbles prevent effective mixing during the MagMAX™ Express-96 run). If bubbles form, quick-spin the plate at 560 × g (~2000 rpm) in a plate centrifuge.
Poor recovery and/or efficiency during spike experiments	<ul style="list-style-type: none"> Incorrectly designed spike amount. <i>or</i> <ul style="list-style-type: none"> Sample concentration is higher than expected. 	Use a spike amount 50–100% of the concentration in the unspiked sample. See Appendix B, “Design guidelines for ProteinSEQ™ System Protein A spike experiments”.
Low spike efficiency	Salt concentration in sample well is too high.	Pre-dilute the sample so that final concentration of salt in the reaction well is <50 mM.
	Matrix interference from IgG or other components.	Evaluate the assay performance with higher sample dilutions.
Increased percent CV	Incorrect plate type used.	Use MagMAX™ Express-96 Skirted Low Profile Plates (Cat. no. 4472783) for Capture, Probes, and Wash plates. See “Label plates” on page 14.



Design guidelines for ProteinSEQ™ System Protein A spike experiments

About spike experiments

Spike recovery is an essential tool for evaluating the accuracy of a quantitation assay in relevant matrices.

A basic spike recovery experiment includes these steps.

1. If necessary, dilute the sample according to the experimental goal.
2. Split the starting sample matrix into two aliquots, one for spiking and a second for referencing.
3. Add a known amount of Protein A (e.g. stock from the standard curve dilution series) into the spike sample, and add a volume of sample diluent (e.g. buffer) equal to the spike volume to the reference sample.
4. Analyze the spike and the reference sample using the same method to generate a mean observed quantitation value.
5. In data processing, subtract the Mean Quantity (Reference) from the Mean Quantity (Spike) to calculate the Reference Adjusted Quantity (RAQ).
6. Divide the RAQ by the Spike Input and multiply by 100 to arrive at a Percent Recovery.




Figure 4 Use of a spike and reference sample to determine percent efficiency for a quantitation assay

Important experimental design considerations

To obtain informative and valid results, consider the following when designing the experiment:

- What concentration of analyte in the spike should be evaluated?
- At what ppm (drug product concentration) should the analyte be evaluated?
- What matrices are available for evaluation?

Guidelines for spike input concentration

The spike input concentration is defined by the goals of the experiment, the concentration of the analyte in the reference sample and the position of that concentration within the standard curve.

The concentration of the spike must be large enough to be differentiated from the analyte concentration already present in the reference sample. Typical spike concentrations range from 50–100% of the reference concentration. Therefore, the reference sample analyte concentrations must be known in order to select the proper spike input concentration.

- If the reference analyte concentration lies between the LLOQ and the mid-point of the standard curve, a 100% spike is recommended.
- If the reference analyte concentration is above the midpoint of the standard curve, a 50% spike is recommended. Note that working at the upper end of standard curve requires care that the final concentration after spiking does not exceed the ULOQ.

Guidelines for matrix selection

The choice of matrix directly affects the design of a spike experiment due to the fact that the ratio of analyte to drug (expressed as ppm in ng analyte/mg drug substance) is a fixed ratio that does not change with dilution. Therefore, practical limitations exist for working with all matrices and a priority must typically be assigned to evaluate either a specific analyte concentration (and the drug substance concentration that follows) or a specific drug substance concentration (and the analyte concentration that follows). For this reason, it is recommended that the goal of the experiment be established followed by procurement of a matrix rather than vice versa.

Volume of sample		250	250	250
Volume of spike		5	10	15
Stock concentrations for spiking	400,000	7,843.14	15,384.62	22,641.50
	40,000	784.31	1,538.46	2,264.15
	8,000	156.86	307.69	452.83
	1,600	31.37	61.54	90.56
	320	6.27	12.30	18.11
	64	1.25	2.46	3.62
	12.8	0.25	0.49	0.72
	2.56	0.05	0.10	0.14
	0.512	0.010	0.019	0.029

Figure 5 Final Protein A concentration of spike using various standard curve stock concentrations and volumes. Green cells = concentrations recommended for spiking studies. Yellow cells = concentrations within the dynamic range but not recommended for spiking studies. Red cells = concentrations out of the ProteinSEQ™ standard curve dynamic range.

Protein A quantitation example

Experimental goals: Evaluate Protein A quantitation at ~150 pg/mL in the presence of ~5 mg/mL IgG. An evaluation matrix containing an estimated 800 pg/mL of residual Protein A and 100 mg/mL IgG is procured.

1. Dilute 5X by mixing 60 μ L sample with 240 μ L sample diluent to reach a volume of 300 μ L.
2. Dilute an additional 2X by adding 150 μ L of 2X Sample Prep reagent to reach a final volume of 600 μ L. The expected concentration is 160 pg/mL Protein A and 20 mg/mL IgG.
3. Generate a pre-reference aliquot of 250 μ L and a pre-spike aliquot of 250 μ L in microcentrifuge tubes, enabling subsequent sample preparation steps. The remaining volume may be discarded.
4. Spike the sample. Figure 5 indicates that the closest spike to 100% can be achieved by spiking 250 μ L of sample with 15 μ L of the 1600 pg/mL stock solution. In order to have a matching reference, 250 μ L of matrix is also combined with 15 μ L of sample diluent to prepare the reference sample. Proceed with the described protocol by boiling, cooling and centrifuging the sample followed by complete transfer of the supernatant to the dilution plate.
5. Run the assay to obtain values for Mean Quantity (Spiked) and Mean Quantity (Reference). In this example, the observed values for Mean Quantity (Spiked) and Mean Quantity (Reference) are 160 pg/mL and 73.6 pg/mL respectively.
6. Calculate the Mean Quantity (Adjusted) for the spiked sample by subtracting Mean Quantity (Reference) from Mean Quantity (Spiked).
7. Calculate Percent Recovery by dividing the Mean Quantity (Adjusted) by the Spike Input and multiplying by 100.




Figure 6 Example Protein A spike experiment with calculations.



Good laboratory practices for PCR and RT-PCR

When preparing samples for PCR or RT-PCR amplification:

- Wear clean gloves and a clean lab coat (not previously worn while handling amplified products or during sample preparation).
- Change gloves whenever you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation and reaction setup.
 - Amplification and analysis of products.
- Do not bring amplified products into the reaction setup area.
- Open and close all sample tubes carefully. Avoid splashing or spraying samples.
- Keep reactions and components capped as much as possible.
- Use a positive-displacement pipettor or aerosol-resistant barrier pipette tips.
- Clean lab benches and equipment periodically with 10% bleach solution or DNAZap™ Solutions (Cat. no. AM9890).



Safety



WARNING! GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
 - Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, etc). To obtain SDSs, see the “Documentation and Support” section in this document.
-

Chemical safety



WARNING! GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the “Documentation and Support” section in this document.
 - Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
 - Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
 - Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended in the SDS.
 - Handle chemical wastes in a fume hood.
 - Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
 - After emptying a waste container, seal it with the cap provided.
 - Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
 - Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
 - **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.
-

Biological hazard safety



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Safety equipment also may include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:
www.cdc.gov/biosafety/publications/bmbl5/BMBL.pdf
- World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:
www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf

Documentation and support

Related documentation

Portable document format (PDF) versions of this guide and the following related documents are available from **thermofisher.com/support**:

Document	Publication number	Description
ProteinSEQ™ Protein A Quantification Kit Quick Reference — Workflow for FAST PCR plates	MAN0013525	Provides information on preparing and running assays using FAST PCR plates.
ProteinSEQ™ Protein A Quantification Kit Quick Reference — Workflow for Standard (non-FAST) PCR plates	MAN0013526	Provides information on preparing and running assays using Standard (non-FAST) PCR plates.

Note: To open the user documentation, use the Adobe™ Reader™ software available from **www.adobe.com**

Note: For additional documentation, see “Customer and technical support” on page 40.

Customer and technical support

Visit **thermofisher.com/support** for the latest in services and support, including:

- Worldwide contact telephone numbers
- Product support, including:
 - Product FAQs
 - Software, patches, and updates
- Order and web support
- Product documentation, including:
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.thermofisher.com/us/en/home/global/terms-and-conditions.html**. If you have any questions, please contact Life Technologies at **www.thermofisher.com/support**.

