

resDNASEQ™ Quantitative DNA Kits

Catalog Numbers 4402085, 4458435, A46014, A26366, A41797, 4464336, 4458441, 4464335 (For CHO, *E. coli*, HEK293, Human, Vero, *Pichia*, NS0, and MDCK)

Pub. No. 4460661 **Rev.** E

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](https://www.thermofisher.com/support).

Product description

The resDNASEQ™ Quantitative DNA Kits are used to quantitate host-cell line residual DNA from CHO, *E. coli*, Human, HEK293, Vero, *Pichia*, NS0, and MDCK cell lines, which are used for production of biopharmaceutical products.

Contents

Table 1 resDNASEQ™ Quantitative CHO DNA Kit (Cat. No. 4402085)

Contents	Amount	Storage
resDNASEQ™ CHO DNA Control		
CHO DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C 2-8°C after first use
resDNASEQ™ CHO Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X CHO DNA Real-Time PCR Assay Mix	300 μL	-25°C to -15°C, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use

Table 2 resDNASEQ™ Quantitative *E. coli* DNA Kit (Cat. No. 4458435)

Contents	Amount	Storage
Box 1, resDNASEQ™ Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X <i>E. coli</i> DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Box 2, resDNASEQ™ <i>E. coli</i> DNA Control		
<i>E. coli</i> DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C before first use 2-8°C after first use

Table 3 resDNASEQ™ Quantitative HEK293 DNA Kit (Cat. No. A46014)

Contents	Amount	Storage
resDNASEQ™ HEK293 DNA Control		
HEK293 DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25 to -15°C before first use 2-8°C after first use
resDNASEQ™ HEK293 DNA Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X HEK293 Assay Mix	300 μL	-25°C to -15°C, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use

Table 4 resDNASEQ™ Quantitative Human DNA Kit (Cat. No. A26366)

Contents	Amount	Storage
resDNASEQ™ Human DNA Control		
Human DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25 to -15°C before first use 2-8°C after first use
resDNASEQ™ Human Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X Human DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use

Table 5 resDNASEQ™ Quantitative MDCK DNA Kit (Cat. No. 4464335)

Contents	Amount	Storage
Box 1, resDNASEQ™ Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X MDCK DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Box 2, resDNASEQ™ MDCK DNA Control		
MDCK DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C before first use 2-8°C after first use

Table 6 resDNASEQ™ Quantitative NS0 DNA Kit (Cat. No. 4458441)

Contents	Amount	Storage
Box 1, resDNASEQ™ Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X NS0 DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Box 2, resDNASEQ™ NS0 DNA Control		
NS0 DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C before first use 2-8°C after first use

Table 7 resDNASEQ™ Quantitative *Pichia* DNA Kit (Cat. No. 4464336)

Contents	Amount	Storage
Box 1, resDNASEQ™ Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X <i>Pichia</i> DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Box 2, resDNASEQ™ <i>Pichia</i> DNA Control		
<i>Pichia</i> DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C before first use 2-8°C after first use

Table 8 resDNASEQ™ Quantitative Vero DNA Kit (Cat. No. A41797)

Contents	Amount	Storage
Box 1, resDNASEQ™ Real-Time PCR Reagents		
2X Environmental Master Mix	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X Vero DNA Assay Mix	300 μL	-25°C to -15°C, protect from light
Box 2, resDNASEQ™ Vero DNA Control		
Vero DNA Control, 30 ng/μL	40 μL	-25°C to -15°C
DNA Dilution Buffer (DDB)	7 mL	-25°C to -15°C before first use 2-8°C after first use

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Life Technologies Ltd | 7 Kingsland Grange | Woolston, Warrington WA1 4SR | United Kingdom
For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4460661

Revision	Date	Description
E	4 November 2019	Add resDNASEQ™ Quantitative HEK293 DNA Kit (Cat. No. A46014).
D	6 February 2019	Update template and legal information.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.