

PrepSEQ[™] Residual DNA Sample Preparation Kit

USER GUIDE

Automated and manual protocols for preparation of samples for
use with resDNASEQ[™] Quantitative DNA Kits

Catalog Number 4413686

Publication Number 4469838

Revision E



Life Technologies Ltd | 7 Kingsland Grange | Woolston, Warrington WA1 4SR | United Kingdom

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4469838

Revision	Date	Description
E	12 Dec 2019	Added the resDNASEQ™ Quantitative HEK293 DNA Kit (Cat. No. A46014).
D	08 June 2018	Updated template, legal, and contact information. Reorganized content. Increase vortex time in elution step to 20 seconds.
C	Dec 2014	Add the resDNASEQ™ Quantitative Human DNA Kit (Cat. No. A26366)

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

TRADEMARKS: All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2019 Thermo Fisher Scientific Inc. All rights reserved.

Contents

■	CHAPTER 1	Product information	5
		Product description	5
		Contents and storage	5
		Required materials not supplied for manual protocols	6
		Required materials not supplied for automated protocols	7
		Workflow	8
■	CHAPTER 2	Prepare the reagents and samples	9
		Prepare the reagents: before first use of the kit	9
		Magnetic beads	9
		Binding Solution	9
		Wash Buffer Concentrate	9
		Prepare reagents: before each use of the kit	9
		Proteinase K (PK) mix	9
		Lysis solution	10
		Guidelines for optimal DNA yields	10
		Sample preparation guidelines	10
		Sample dilution (if necessary)	10
		Triplicate extractions	11
		Extraction control guidelines	11
■	CHAPTER 3	Manual protocol for residual DNA extraction	13
		Digest the test samples and controls	13
		Bind the DNA	13
		Wash the DNA	14
		Elute the DNA	15
■	CHAPTER 4	Automated protocol for residual DNA extraction	16
		Before each use of the kit	16
		Ensure that you have the correct plates	16
		Prepare the plates	16
		Prepare the lysis plate	16
		Process samples on the instrument	17

- **APPENDIX A** Troubleshooting 19
- **APPENDIX B** Good laboratory practices 20
 - Good laboratory practices for PCR and RT-PCR 20
 - Avoiding false positives due to cross-contamination 20
 - Plate layout suggestions 20
- **APPENDIX C** Safety 21
 - Chemical safety 22
 - Biological hazard safety 23
- Documentation and support 24**
 - Related documentation 24
 - Customer and technical support 24
 - Limited product warranty 25



Product information

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Product description

The PrepSEQ™ Residual DNA Sample Preparation Kit extracts host-cell DNA from products that are produced in cell lines such as CHO, *E. coli*, HEK293, Human, Vero, *Pichia*, NS0, and MDCK. The kit uses chemical lysis and magnetic beads to extract genomic DNA from diverse sample types, including samples that contain high protein and low DNA concentration.

For quantification of host-cell line residual DNA, we recommend using the resDNASEQ™ Quantitative DNA Kits as described in the *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836). To ensure accurate quantitative results, extract each sample in triplicate and perform a single PCR reaction for each extraction.

Contents and storage

The kit contains reagents sufficient for 100 extractions.

Table 1 PrepSEQ™ Residual DNA Sample Preparation Kit (Cat. No. 4413686)

Contents	Amount	Storage
Box 1, PrepSEQ™ Nucleic Acid Extraction Kit		
Lysis Buffer	2 × 50 mL	Room temperature
Binding Solution (Isopropanol), empty bottle	1	
Wash Buffer Concentrate	2 × 26 mL	
Elution Buffer	25 mL	
Proteinase K (PK) Buffer Can be used for existing validated manual protocols.	50 mL	
Proteinase K (PK) Buffer II ^[1] Recommended for new manual protocols. Required for automated protocols.	50 mL	

Contents	Amount	Storage
Box 2, PrepSEQ™ Nucleic Acid Extraction Kit		
Magnetic Particles	2 × 1.5 mL	Room temperature
Box 3, PrepSEQ™ Nucleic Acid Extraction Kit		
Proteinase K, 20 mg/mL	1.25 mL	–20°C or below
PrepSEQ™ Residual DNA Sample Preparation Kit		
Proteinase K, 20 mg/mL	1.25 mL	–20°C or below
Yeast tRNA, 10 mg/mL	0.5 mL	
Glycogen, 5 mg/mL	2 × 1.0 mL	

^[1] Also sold separately (Cat. No. 4415320).

Required materials not supplied for manual protocols

Unless otherwise indicated, all materials are available through **thermofisher.com**.
 MLS: Fisher Scientific (**fisherscientific.com**) or other major laboratory supplier.

Item	Source
Equipment	
Pharma Magnetic Stand-96 or Magnetic stand, 16-position	A31543 12321D
Block heater for use with 2-mL tubes. Manual DNA extraction involves two incubations at different settings, so two heaters may be convenient.	Major laboratory supplier (MLS)
Benchtop microcentrifuge for 1.5-mL and 2-mL tubes	MLS
Vortex-Genie™ 2T Mixer	VWR™ 14216-188, VWR™ 14216-186
Vortex Adapter-60, for use with the Vortex-Genie™	AM10014
Consumables	
Disposable gloves	MLS
Aerosol-resistant micropipette tips	MLS
PIPETMAN™ Pipettors, P1000, P200, P20 and P10: <ul style="list-style-type: none"> • Positive-displacement • Air-displacement • Multichannel 	MLS

Item	Source
Nonstick, RNaseZap™-free Microfuge Tubes, 1.5-mL (1 box; 250 tubes/box)	AM12450
Safe-Lock PCR clean microcentrifuge tubes, round-bottom, 2-mL	VWR™ 62111-754
Reagents	
Ethanol, 95% IMPORTANT! Do not use denatured ethanol. It contains components that are not compatible with the protocol.	MLS
Isopropanol, 100%	MLS
5 M NaCl and 1 N NaOH solutions	MLS
Hydrochloric acid (HCl)	MLS
PBS solution, 1X, pH 7.4, free of Mg and Ca (if needed to dilute samples before DNA extraction)	14190094

Required materials not supplied for automated protocols

Unless otherwise indicated, all materials are available through **thermofisher.com**.
MLS: Fisher Scientific (**fisherscientific.com**) or other major laboratory supplier.

Table 2 Pharma KingFisher™ Flex-96 Deep Well Magnetic Particle Processor (Cat. No. A31508) and Pharma MagMAX™ Express-96 instrument^[1] accessories

Item	Source
Pharma MagMAX™ Express-96 DW plate	A31540
Pharma MagMAX™ Express-96 Deep-Well Tip Combs	A31537
Pharma KingFisher™ Flex Magnetic Head for 96 Deep-Well Plate	A31542
Pharma MagMAX™ 96 PCR Well Magnetic Head	4472991
Pharma MagMAX™ Express-96 Standard Plates	A31541

^[1] This instrument is no longer available for purchase

Table 3 Additional materials

Item	Source
Consumables	
Disposable gloves	MLS
Aerosol-resistant micropipette tips	MLS

Item	Source
PIPETMAN™ Pipettors, P1000, P200, P20 and P10: <ul style="list-style-type: none"> • Positive-displacement • Air-displacement • Multichannel 	MLS
Nonstick, RNaseZap™-free Microfuge Tubes, 1.5-mL (1 box; 250 tubes/box)	AM12450
Reagents	
Ethanol, 95% IMPORTANT! Do not use denatured ethanol. It contains components that are not compatible with the protocol.	MLS
Isopropanol, 100%	MLS
5 M NaCl and 1 N NaOH solutions	MLS
Hydrochloric acid (HCl)	MLS
PBS solution, 1X, pH 7.4, free of Mg and Ca (if needed to dilute samples before DNA extraction)	14190094

Workflow

Manual residual DNA extraction	Automated residual DNA extraction
Digest the test samples and controls (page 13)	Prepare the plates (page 16)
▼	▼
Bind the DNA (page 13)	Prepare the lysis plate (page 16)
▼	▼
Wash the DNA (page 14)	Process samples on the instrument (page 17)
▼	
Elute the DNA (page 15)	

2

Prepare the reagents and samples

Prepare the reagents: before first use of the kit

Magnetic beads

1. Set a block heater to 37°C.
2. Incubate the Magnetic Particle suspension at 37°C for a minimum of 10 minutes with intermittent vortexing at setting #7, or until the particles are completely suspended.

Binding Solution

1. Add 30 mL of 100% isopropanol to the Binding Solution bottle.
2. Label the bottle to indicate that it contains isopropanol, then store the bottle at ambient temperature.

Wash Buffer Concentrate

1. Add 74 mL of 95% ethanol to one bottle of PrepSEQ™ Wash Buffer Concentrate, then mix completely.
2. Label the bottle to indicate that it contains ethanol, then store the bottle at room temperature.

Prepare reagents: before each use of the kit

Proteinase K (PK) mix

- Use Proteinase K (PK) Buffer II for new manual protocols and automated protocols.
- Note:** Proteinase K (PK) Buffer is also provided in the kit for use with existing manual protocols that have been internally validated with this buffer.
- Prepare a fresh mix before each use of the kit.
- Include a 10% overage to account for pipetting losses.

Component	Number of extractions				
	1	7	10	13	25
Proteinase K, 20 mg/mL	10 µL	70 µL	100 µL	130 µL	250 µL
Proteinase K (PK) Buffer II or Proteinase K (PK) Buffer	60 µL	420 µL	600 µL	780 µL	1,500 µL

Lysis solution

- Prepare a fresh mixture immediately before use or during Proteinase K incubation.
- Prepare 360 μL (amount required) of lysis solution mix per sample.

Reagent	Volume for ~20 extractions
Glycogen, 5 mg/mL	180 μL
Yeast tRNA, 10 mg/mL	4 μL
Lysis Buffer	7,600 μL
Total	7,784 μL

Guidelines for optimal DNA yields

- Maintain a homogenous suspension of the magnetic beads to maximize the surface area to which the DNA can bind. The appearance of the mixture should be homogenous after mixing.
- After drying, the DNA remains bound to the magnetic beads. Do not allow the magnetic beads to over-dry because this reduces the elution efficiency; over-dried beads are not easily resuspended.
- During manual elution, vortex every 2 minutes to assist elution. This will result in better yield during recovery.

Note: Some test samples cause the beads to adhere very firmly to the tube wall, while others form loose pellets that detach during the vortex steps. All pellets should dissolve with vortexing during heated elution. If vortexing does not result in full resuspension, then wash the beads off the tube by pipetting.

Note: White or brown precipitate may form in the Magnetic Particles tube if it is stored at 2–8°C. The precipitate will dissolve when it is heated to 37°C for a minimum of 10 minutes with intermittent vortexing. Make sure the precipitate is completely dissolved before using the beads.

Sample preparation guidelines

Sample dilution (if necessary)

Test samples from the early purification process often contain levels of DNA that are above the highest point of the residual DNA assay standard curve. You must dilute these samples (from 1:100 up to 1:10,000) before PrepSEQ™ Residual DNA sample preparation.

- Dilute test samples before DNA extraction with a solution of 1X PBS (pH 7.4; free of Mg and Ca) or 50 mM Tris, pH 8.0, 0.5 M NaCl.

Note: Diluting samples in water or TE reduces extraction efficiency.

- Use the sample dilution buffer as the negative extraction control instead of water.
- Alternatively, dilute extracted DNA with elution buffer before running the PCR reaction.

Triplicate extractions

Triplicate extractions are required for post-PCR analysis calculation of mean quantity, standard deviation, and coefficient of variation.

In addition to test samples, we recommend triplicate extractions for controls (for an explanation of controls, see “Extraction control guidelines” on page 11).

Perform a single PCR reaction for each extraction.

The table below illustrates the total number of extractions required based on the 1, 2, and 3 samples extracted in a batch.

Table 4 Total number of extractions per batch of test samples

Number of test samples		Total number of extractions for the batch
1	3 extractions required for each: <ul style="list-style-type: none"> • Test sample • Test sample extraction/recovery control (ERC) • Negative extraction control^[1] 	9
2		15
3		21

^[1] Optional during routine testing.

Extraction control guidelines

We recommend that you use the following extraction controls:

Type of control	Contains	Number to run	Used to
Negative (NEG) ^[1]	1X PBS	1 per batch of extractions	Test for cross-contamination of DNA extraction reagents.
Extraction/recovery (ERC)	Positive control from the resDNASEQ™ Quantitative DNA Kit	3 per sample	<ul style="list-style-type: none"> • Evaluate the efficiency of DNA extraction, recovery, and quantification from test samples. • Verify reagent and system performance.

^[1] Optional during routine testing.

For the Extraction/recovery (ERC) :

- Prepare the control standard dilutions as described in *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836).
- Add a volume of positive control standard dilution to each test sample that yields a PCR input DNA amount that is 2–10 times the amount of DNA measured in the test sample *without* the addition of the DNA control.

For example:

- The DNA amount measured in a test sample is ≤ 1 pg.
- To prepare a 10 pg ERC for a PCR elution volume of 50 μL , spike samples with 16.7 μL of the 3 pg/ μL positive control standard dilution (SD3) = 50 pg spike to yield 10 pg per PCR reaction.
- Prepare three separate extractions for each test sample, then add the ERC to each reaction. Do not prepare a large volume of ERC, then aliquot it into three reactions.

Note: To calculate the efficiency of DNA recovery and quantification from the test samples, subtract the amount of DNA measured in the sample *without* the addition of DNA control from the amount of DNA measured in the ERC sample.



Manual protocol for residual DNA extraction

Digest the test samples and controls

1. Set a block heater to 56°C. If available, set a second block heater to 70°C.
2. Label 2-mL Safe-Lock tubes:
 - 3 for each sample
 - 3 for each sample + ERC
 - 3 for NEG
3. Add 100 µL of sample, sample + ERC, or 1X PBS to into each tube.
4. Add 10 µL of 5 M NaCl and 70 µL Proteinase K/Proteinase K Buffer II mix.
5. Briefly vortex and centrifuge.
6. Incubate at 56°C for 30 minutes.

If only one block heater is available, after this incubation step is complete, reset the block heater to 70°C for the elution step.

Note: For samples with high protein concentration, extending the incubation time to 60 minutes can increase recovery.
7. Cool samples to room temperature.
8. Add 360 µL freshly made Lysis solution mix to each tube.

Bind the DNA

1. Vortex the Magnetic Particles to resuspend the particles.

Note: The appearance of the mixture should be homogeneous.
2. Add 30 µL of the Magnetic Particles to each tube.
3. Add 400 µL Binding Solution to each tube.
4. Mix and vortex the tubes:
 - a. Close the caps, **immediately** invert each tube twice to mix.

- b. Vortex the tubes in the vortex adaptor for 5 minutes at setting #7.
5. Briefly centrifuge the tubes for 15 seconds at top speed ($>15,000 \times g$) to collect the Magnetic Particles at the bottom of the tubes.
6. Place the tubes in the magnetic stand with the pellet against the magnet, then let the tubes stand for 5 minutes or until the solution is clear.
7. Without disturbing the magnetic beads, remove the supernatant using a PIPETMAN™ pipette or by aspiration.

Wash the DNA

For aspiration of liquid supernatants and wash buffers during sample preparation, we recommend attaching the waste-collection bottle to the vacuum using tubing that can accommodate 200- μ L pipette tips.

1. Remove the tube rack (with tubes) from the magnetic stand, then add 300 μ L of Wash Solution to the tubes. Vortex the tubes for 5 seconds at room temperature at setting #7.
2. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
3. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA are magnetically captured after approximately 1 minute.
4. Without disturbing the Magnetic Particles, remove the supernatant using a PIPETMAN™ pipette or by aspiration.
5. Remove the tube rack (with tubes) from the magnetic stand, then add 300 μ L of Wash Solution to each tube for a second wash. Vortex the tubes for 5 seconds at room temperature at setting #7.
6. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
7. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA are magnetically captured after approximately 1 minute.
8. Open all tubes, then start the 5-minute timer.
9. Without disturbing the Magnetic Particles, remove the supernatant using a PIPETMAN™ pipette or by aspiration.
Use a P200 to remove the remaining solution from the bottom of the tube.

10. With the tube lid open, air-dry the Magnetic Particles pellet in the magnetic stand for no more than 5 minutes at room temperature.

IMPORTANT! Air-dry to remove ethanol from the Wash Solution. After dried, the DNA stays bound to the magnetic beads. Do not over-dry; over-dried beads are not easily resuspended.

Elute the DNA

1. Add 50 μ L of Elution Buffer to each tube.
2. Vortex the tubes for 20 seconds at high speed, then incubate the tubes at 70°C for 7 minutes. Vortex the tubes two to three times during the incubation to help resuspension.

Note: (Optional) If vortexing does not result in full resuspension, then wash the beads off the tube by pipetting.

3. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
4. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA are magnetically captured after approximately 1 minute.
5. Without disturbing the Magnetic Particles, transfer the liquid phase containing the eluted DNA to a new nonstick 1.5-mL microcentrifuge tube.
6. Centrifuge the tube at top speed ($>15,000 \times g$) for 3 minutes to collect the Magnetic Particles at the bottom, then place the tubes in the magnetic stand for 1 minute.
7. Without disturbing the Magnetic Particles, transfer the liquid phase containing the eluted DNA to the nonstick 1.5-mL microcentrifuge tube.

Note: Magnetic Particles can inhibit PCR.

When you finish the sample extraction procedure, see the *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836) to set up PCR for DNA quantification.

Use 10 μ L of the eluate in the real-time PCR.



Automated protocol for residual DNA extraction

You can use the KingFisher™ Flex or MagMAX™ Express 96-deep well automation platforms to automate the extraction of host-cell line residual DNA. For all chemicals, read the Safety Data Sheet (SDS) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Before each use of the kit

Ensure that you have the correct plates

The KingFisher™ Flex or the MagMAX™ Express require 5 plates.

Plate name	Plate type
Lysis	96 deep-well plate
Wash 1	96 deep-well plate
Wash 2	96 deep-well plate
Elution	96 deep-well plate
Comb loading plate	96 deep-well tip comb combined with 96 standard plate

Prepare the plates

Prepare the Wash 1, Wash 2, and Elution plates:

Plate name	Plate type	Volume of buffer to add
Wash 1	96 deep-well plate	300 µL of Wash buffer
Wash 2	96 deep-well plate	300 µL of Wash buffer
Elution	96 deep-well plate	200 µL of Elution buffer

Prepare the lysis plate

In all steps that require pipetting, dispense liquid at bottom center of the wells.

1. Add 100 µL to the appropriate wells of the 96 deep-well Lysis plate:
 - 3 wells for each sample
 - 3 wells for each sample + ERC
 - 3 wells for NEG

2. Add 10 µL of 5 M NaCl to each sample well.
3. Add 70 µL Proteinase K/Proteinase K (PK) Buffer II mix to each sample well.

Process samples on the instrument

1. Select the script or program for the instrument you are using:

Instrument	Select
KingFisher™ Flex	PrepSEQ_resDNA_v1 script
MagMAX™ Express-96	PrepSEQ_resDNA_2011 PrepSEQ_1hr_resDNA (if installed)

2. Load the plates into the instrument in the order listed below. After loading each plate, press **START** to move the turntable.
 - a. Comb loading plate
 - b. Elution plate with 200 µL of Elution Buffer
 - c. Wash 2 plate with 300 µL of wash buffer
 - d. Wash 1 plate with 300 µL of wash buffer
 - e. Lysis plate
3. Press **START** to begin the PK digestion process.
The instrument mixes the samples for 10 seconds at fast speed, then incubates the samples at 56°C for 30 minutes, mixing at slow speed. When digestion is complete, the instrument pauses and returns the Lysis plate to the loading position.
4. After the digestion step is complete, add additional components to the Lysis plate:
 - a. Remove the Lysis plate from the instrument.
 - b. Add 360 µL of Lysis Solution to each sample well.
 - c. Add 30 µL of Magnetic Particle suspension to each sample well.
 - d. Add 400 µL of Binding Solution to each sample well, then immediately pipet up-and-down three times to mix.
 - e. Place the plate back into the instrument loading position, then press **START** to begin binding.
5. When DNA extraction is finished, the instrument returns the Elution plate to the loading position.

When you finish the sample extraction procedure, refer to the *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836) to set up PCR for DNA quantitation.

Note: Use 10 μ L of the eluate in the real-time PCR.



Troubleshooting

Observation	Possible cause	Action
Poor extraction efficiency (low yields)	Overdrying the sample.	Start the 5-minute timer before removing ~300 μ L from the first 6–8 samples. Then continue removing wash buffer from the remaining samples.
	Magnetic Particles are difficult to resuspend during the elution.	Incubate the pellets at 70°C for > 7 minutes. Vigorously vortex the tubes three times during incubation to help resuspension. Do not overdry. If necessary, repeat the incubation and vortexing steps.
	Formation of precipitate in Magnetic Particles.	Incubate the Magnetic Particle suspension at 37°C with intermittent vortexing at setting #7 until the particles are completely suspended.
	PK Buffer was used instead of PKII Buffer.	Use PKII Buffer.
Particles no longer produce consistent results (fine brown sandy particles and brown color are observed in the supernatant)	Samples have low pH.	Measure the pH of the sample and adjust the pH to between 6 and 8.
	Magnetic Particles were stored at –20°C.	Order new materials and store them at room temperature.



Good laboratory practices

Good laboratory practices for PCR and RT-PCR

- Wear clean gloves and a clean lab coat.
 - Do not wear the same gloves and lab coat that you have previously used when handling amplified products or preparing samples.
- Change gloves if you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation and reaction setup.
 - Amplification and analysis of products.
- Do not bring amplified products into the reaction setup area.
- Open and close all sample tubes carefully. Avoid splashing or spraying samples.
- Keep reactions and components capped as much as possible.
- Use a positive-displacement pipettor or aerosol-resistant barrier pipette tips.
- Clean lab benches and equipment periodically with 10% bleach solution or DNA decontamination solution.

Avoiding false positives due to cross-contamination

To avoid false positives due to cross-contamination:

- Prepare and close all negative control and unknown sample tubes before pipetting the positive control.
- Do not open tubes after amplification.
- Use different sets of pipettors when pipetting negative control, unknown, and positive control samples.

Plate layout suggestions

- For each plate row, dispense in sequence from left to right the: negative controls, unknown samples and ERCs, and positive controls (at the end of the row or column).
- Place positive controls in one of the outer columns (10–12).
- If possible, separate all samples from each other by at least one well; if space is limiting, place at least one well between unknown samples and controls.



Safety



WARNING! GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
 - Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, see the “Documentation and Support” section in this document.
-

Chemical safety



WARNING! GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.



WARNING! HAZARDOUS WASTE (from instruments). Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.



WARNING! 4L Reagent and Waste Bottle Safety. Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.



Biological hazard safety



WARNING! Potential Biohazard. Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:
<https://www.cdc.gov/labs/pdf/CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2009-P.pdf>
 - World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:
www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf
-

Documentation and support

Related documentation

Document	Publication number	Description
<i>PrepSEQ™ Residual DNA Sample Preparation Kit Quick Reference</i>	4469839	For brief instructions on using the PrepSEQ™ Residual DNA Sample Preparation Kit.
<i>resDNASEQ™ Quantitative DNA Kits User Guide</i>	4469836	For information on performing PCR after sample extraction.
<i>resDNASEQ™ Quantitative DNA Kits Quick Reference</i>	4469837	For brief instructions on using the resDNASEQ™ Quantitative DNA Kits.
<i>Thermo Scientific™ KingFisher™ Flex User Manual</i>	N07669	For information on the KingFisher™ Flex instrument.
<i>Applied Biosystems™ MagMAX™ Express 96 User Manual</i>	N07848	For information on the MagMAX™ Express 96 DW instrument.

Customer and technical support

Visit thermofisher.com/support for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments
- Order and web support
- Product documentation
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

