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1. 使用细胞工厂CF-40, 无血清瞬时表达抗体
Serum-free Transient Expression of Antibodies Using 40 Tray Cell Factories
Peter Schlenke et al, Poster Renschler Biotechnologie
HEK 293 cells, antibody expression, serum free
2. 使用NUNC细胞工厂, T-flasks和培养转瓶(配合HyQ SFM4Megavir), Vero细胞培养的脊髓灰质炎病毒产品
Production of Poliovirus using Vero Cells Cultured in NUNC Cell Factories, T-flasks and Roller Bottles with HyQ SFM4Megavir. HyClone Application Note
Polio virus, VERO cells, SFM4MegaVir
3. 使用NUNC细胞工厂在封闭系统中生产大量树突状细胞
Generation of large number of dendritic cells in a closed system using Cell Factories. Sandra Tuyaerts et al, JIM 264 (2002) 135-151
Dendritic cells
4. 使用NUNC细胞工厂的腺病毒大规模制备
Large scale preparation of adenovirus using NUNC Cell Factories. work protocol
Adenovirus, 293 cells
5. 使用NUNC细胞工厂, 成熟的单核细胞来源树突状细胞(应用于临床)的大规模制备
Large scale generation of mature monocyte-derived dendritic cells for clinical application in Cell Factories. Thomas Berger et al, JIM 268 (2002) 131-140
Dendritic cells
6. Lentiviral基因转入初级和二级NOD/SCID再植细胞
Lentiviral gene transfer into primary and secondary NOD/SCID repopulating cells. Niels-Bjarne Woods et al. Blood, Vol 96, No.12, pp. 3725-3733
293T cells, Lentivirus
7. VHS病毒(Viral Hemorrhagic Septicemia)分子内二硫键的鉴定和VHS病毒糖蛋白二硫键的鉴别。VHS病毒为一种鱼类弹状病毒
Characterization of Intramolecular Disulfid Bonds and secondary Modifications of the glycoprotein from Viral Hemorrhagic Septicemia Virus, a Fish Rhabdovirus. Katja Einer-Jensen et al, Journal of Virology 1998, p. 10189-10196, Vol. 72, No.12
BF-2 cells, Hemorrhagic Septicemia Virus, Fish Rhabdovirus
8. 通过免疫黑色素瘤患者, 接种成熟的、冻存的、多肽致敏单核细胞来源树突状细胞, 从而快速介入肿瘤特异性型辅助性T淋巴细胞
Rapid Introduction of Tumor specific Type 1T Helper cells in metastatic Melanoma patients by vaccination with mature, cryopreserved, peptide-loaded monocyte-derived dendritic cells. Beatrice Schuler-Thurner et al
Dendritic cells
9. Renschler Biotechnologie, Corporate Presentation, 2004五月(报告)
Renschler Biotechnologie, Corporate Presentation, May 2004
Contract manufacturing
10. 使用Q Sepharose XL, 腺病毒的快速纯化
Rapid Adenovirus Purification using Q Sepharose XL. Kjell Eriksson et al, poster Wilbio, Lake Tahoe Nov 2001 + ESGT, Stockholm Oct.2000
Adeno virus, gene therapy
11. 一种适用于肾小管上皮细胞生长及病毒生产的通用无血清培养基
A versatile serum-free medium for kidney epithelial cell growth and virus production Paul Price et al, Gibco, Focus 2002, vol 24, p. 24-28
MDCK cells, MDBK cells, VERO cells, PK-15 cells, BHK 21 cells
12. 产品标书 (Richard Snyder, 佛罗里达大学)
Production Bidding Form University of Florida, Richard Snyder
293 cells, vector production
13. GMP-MSK的临床规模化生产
GMP – Clinical scale production of MSC L.Sensebe, presentation Tissue Engineering 2004, Bordeaux
MSC, bone marrow cells
14. 成人骨髓基质干细胞体外软骨形成, 定义软骨发生过程中的细胞和分子事件顺序
In vitro cartilage formation by human adult stem cells from bone marrow stroma defines the sequence of cellular and molecular events during chondrogenesis. Ichiro Sekiya et al. PNAS, april 2, 2002, vol 99, no.7, 4397-4402
Human marrow stromal cells, MSC
15. MRI磁共振成像、靶向导管定位移植成肌祖细胞至梗死左心室心肌
Magnetic Resonance Imaging of targeted catheter-based Implantation of myogenic precursor cells into infarcted left ventricular Myocardium. Jerome Garot et al., Journal of the American College of Cardiology, Vol 41, No.10, p. 1841-1846, 2003
Myogenic precursor cells (MPC)
16. 核仁素是一种基质结合区DNA结合蛋白, 能够特异地识别潜在的碱基未配对高发区
Nucleolin is a Matrix attachment region DNA-binding Protein that specifically recognizes a region with high base-unpaired Potential. Lillian Dickinson et al., Molecular and cellular Biology, Jan 15, 1995, vol.15, no.1, p. 456-465
K562 (human erythroleukemia cell line)
17. 人胸腺基质淋巴生成素有效刺激髓细胞
Human Thymic Stromal Lymphopoietin preferentially stimulates Myeloid cells. Pedro Reche et al., Journal of Immunology
Human epithelial cells, 293T
18. 人眼前节内第二代FIV载体介导标记基因的表达式后房水流动性维持
Preservation of aqueous outflow facility after second generation FIV vector mediated expression of marker genes in anterior segments of human eyes. Niels Loewen et al. Investigative Ophthalmology visual science, Dec. 2002, vol. 43, No.12, 3686-3690
Lentiviral vectors, 293T cells
19. 膜辅蛋白是流行性角膜炎关联腺病毒的一个受体
Membrane cofactor protein is a receptor for adenoviruses associated with epidemic eratoconjunctivitis. Eugene Wu et al., Journal of Virology, Apr. 2004, p. 3897-3905
Chang conjunctival cells

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Immune response to Hepatitis A virus capsid proteins after infection. Chwan-Heng Wang et al., Journal of Clinical Microbiology, Mar. 1996, Vol. 34, No.3

Fibroblast cells, Hepatitis A virus

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Myoblasts transplanted into rat infarcted myocardium are functionally isolated from their host. Bertrand Leobon et al., PNAS, June 24, 2003, Vol.100, No.13 7808-7811

Primary muscle cells (Wistar rat), transfection, Adenovirus

22. 凝血因子VII突变体V154G，模拟凝血因子VIIa的一种酶原样形态
Factor VII mutant V154G models a zymogen-like form of Factor VIIa. Raffaella Toso et al., Biochemical Journal Immediate Publication, Oct. 1, 2002

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The N-terminal globular domain of the laminin a1 chain binds to a1b1 and a2b1 integrins and to the heparin sulfate-containing domains of perlecan. Norbert Ettner et al., FEBS Letters 430 (1998) 217-221

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Theoretical comparing culture systems

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Phosphoglycerate kinase acts in tumor angiogenesis as a disulphide reductase. Angelina Lay et al., Nature Vol. 408 Dec. 14, 2000 869-873

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Immortalization and characterization of a cell line exhibiting a severe multiple sulphatase deficiency phenotype. Kathy Nelson et al., Biochem J. (1997) 326 (125-130)

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Rapid expansion of recycling stem cells in cultures of plastic-adherent cells from human bone marrow. David Colter et al., Proc Natl Acad Sci USA 2000 March 28; 97 (7): 3213-3218

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Expression, purification and characterization of recombinant caprine N-acetylglucosamine-6-sulphatase. (Tom Litjens et al. Biochem J. (1997) 327 (89-94))

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In vitro generated cytolytic T lymphocytes reactive against head and neck cancer recognize multiple epitopes presented by HLA-A2, including peptides derived from the p53 and MDM-2 proteins. Tadao Asai et al., Cancer Immunity, Vol 2, p.3 (April 16, 2002)

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Expansion of human adult stem cells from bone marrow stroma: Conditions that maximize the yields of early progenitors and evaluate their quality. Ichito Sekiya et al., Stem Cells 2002; 20: 530-541

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Expression, purification and reconstitution of receptor for pituitary adenylate cyclase-activating Polypeptide. Tetsuya Ohtaki et al., J. Biol Chem, Vol 273, issue 25, 15464-15473 (1998)

PACR19 cells

34. US专利6,664,385/ 代理: Biogen US Patent 6,664,385/ Assignee: Biogen

Fusion protein KIM-1, COS cells

35. 应用组合C4S-Keyhole成像技术来研究冷冻保护剂在工程组织中的动力学
Use of a combined C4S-Keyhole imaging technique to study the dynamics of cryoprotective agents in an engineered tissue. N.P.Bideault et al., Presentation at ISMRM 98, Sydney Australia

Fibroblasts

36. 使用Daniplestim、Leridistim、Progenipoiectin、Promegapoiectin和自体血浆进行粒系祖细胞和粒系祖后细胞的临床规模生产
Clinical scale production of granulocyte progenitor and post-progenitor cells using daniplestim, leridistim, Progenipoiectin, Promegapoiectin and autologous plasma. S.Patel et al., Cytotherapy 2000,vol2, no.2, pp.85-94

PBMC

37. 纤维蛋白原通过Aa-羧基端RGD位点结合整合蛋白a5b1
Fibrinogen binds to Integrin $\alpha 5 \beta 1$ via the carboxyl-terminal RGD site of the Aa-chain. Kazuhisa Suehiro et al., J.Biochem. Vol.128, pp.705-710 (2000)

Recombinant Fibrinogen, transfected BHK cells

38. US专利: 6,723,325/ 代理: Acambis

US Patent: 6,723,325 Assignee: Acambis

Smallpox vaccine, vaccinia virus, MRC-5, chicken embryo fibroblasts

39. 使用NUNC细胞工厂对粒系祖细胞和粒系祖后细胞进行的临床规模生产
Clinical scale production of granulocyte progenitors and post-progenitors using NUNC Cell Factories. R. Guo et al., ISHAGE, San Diego 2000, abstract

Peripheral blood progenitor cells

40. 对用于生产人乙酰胆碱酯酶（使用重组293细胞）的锚着依赖细胞繁殖系统的评估
Evaluation of anchorage-dependent cell propagation systems for production of human acetylcholinesterase by recombinant 293 cells. A. Lazar et al., Cytotechnology, 1993;13(2): 115-123

Recominant 293 cells

41. 在鸟类细胞中表达外源基因：一种潜在的重组蛋白生产工具
Heterologous Gene Expression in Avian Cells: Potential as a producer of recombinant (Acrobate) Proteins. Sun-Young Lee et al., J.Biomed.Sci 1999; 6:8-17

QF cell line, EPO expression

42. 从细胞培养中快速高效纯化甲肝病毒
Rapid and efficient purification of Hepatitis A virus from cell culture. Naomi Bishop et al, JIM 47 (1994) 203-216

Hepatitis A virus, BS-C-1 cells (African green monkey kidney)

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Large scale transfection of mammalian cells for the production of recombinant adenovirus. Shyam Subramanian, Merck&Co, presentation Wilbio November 2000 Gene Therapy

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Optimization of retroviral vector production using serum free conditions. H.G.Eckert et al, Eufets, workpaper in collaboration with NUNC

Serum-free conditions, retroviral vectors, CD34+ cells

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Incubator design for optimal heat transfer and temperature control in plastic bioreactors. William Adams, et al. Merck&Co, Pharmaceutical engineering Jan/feb 2003, Vol 23, no. 1

Heat transfer, Temperature control

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Toxoplasma gondii in vitro cultivation: economic and efficient mass production. C. Harmer et al., Journal of Microbiological Methods 27 (1996), 225-228

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Purification and characterization of glutathione S-transferase from the human hepatoma derived PLC/PRF/5 cell line. Paul Dierickx, Biomedical research 10 (4) 301-306, (1989)

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Large scale Production of Hepatitis A virus in cell culture: Effect of type of infection on virus yield and cell integrity. Betty Robertson et al, J. GEN. Virol. (1988), 69, 2129-2134

Hepatitis A virus, FRhK4 cells, foetal resus monkey kidney cells

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Large scale production of human fibroblast interferon in multitray battery system. V. Pakos et al. Develop.biol.Standard., vol. 60 pp. 317-320

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Communicating vessel systems for mass cell culture of anchorage dependent cells. R. Skoda et al., Develop.biol.Standard. 42, 1978, 121-126

General handling description

51. 甲肝病毒的体外繁殖和检验
Propagation and assay of hepatitis A virus in vitro. Guenter Siegel et al., Journal of Virological Methods, 9 (1984) 53-67
Hepatitis A, PLC/PRF/5 cells (hepatoma derived), MRC-5 cells
52. 用于恶性肿瘤过继免疫治疗的LAK细胞的新型大规模生产和冷冻工艺流程
A new procedure for large scale production and freezing of lymphokine activated killer (LAK) cells to be used in adoptive immunotherapy of cancer. Carlo Gambacorti et al. Tumori, 74: 523-530 (1988)
PBL activation
53. 肯塔基州农村儿童接种的一种福尔马林灭活甲肝病毒疫苗 (VAQTA), 安全、耐受且具免疫原性
Safety, tolerability and immunogenicity of a formalin-inactivated hepatitis A vaccine (VAQTA) in rural Kentucky children. Stan Block et al. Pediatr. Infect. Dis. J. 1993, 12:976-980
MRC-5 cells, Hepatitis A
54. 琼脂糖胶电泳系统分离抗体在动物农业学中的应用
Agarose Gel Electrophoresis System for the separation of Antibiotics used in animal Agriculture. Michael Salvatore et al., Analyst march 1993, vol 118
CF as Bioassay dish
55. 重组小鼠神经营养因子3 (neurotrophin-3) 的生产和鉴定
Production and characterization of recombinant mouse neurotrophin-3. Rudolf Goetz et al., Eur.J.Biochem. 204, 745-749 (1992)
Rabbit kidney cells, RK13, vaccinia virus
56. 人血管生长素 (Angiogenin) 在幼仓鼠肾细胞中的表达
Expression of human Angiogenin in cultured baby hamster kidney cells. Kotoku Kurachi et al. Biochemistry 1988, 27, 6557-6562
BHK cells, transfection
57. 一个核酸酶在一种甲肝病毒疫苗VAQTA纯化过程中的应用
Use of a nuclease enzyme in the purification of VAQTA, A hepatitis A vaccine A. Hagen et al., Biotechnology and applied Biochemistry 23(3), 1996, 209-215
Large scale manufacture, inactivated Hepatitis A vaccine
58. 重组人高血压蛋白原酶的酶活性鉴定
Enzymatic characterization of recombinant human rennin. Pilot et al. Biochemistry and Cell Biology 73 (3-4). 1995. 163-170
Transfected DAMP cells, dog pithelial cells, human rennin
59. 用于细胞培养的非灌注附加系统
Nonperfused attachment systems for cell cultivation. AY Elliot, Bioprocess Technol (United States) 1990, 10, p207-116; ISSN 0888-7470
Large scale application description
60. Mammoduline的大规模生产; 发展一种贴附细胞生长的低成本高产工艺
Large scale production of mammoduline. Development of a low cost, large scale production process for adherent cell growth. EJ Schlaeger et al., ESACT-Proceedings 1994
MDA-MB-21 cells, CF as carrier-seeding device
61. 在特定蛋白缺失培养基中培养的人癌细胞肿瘤分泌产物的分离
Isolation of tumor-secreted products from human carcinoma cells maintained in a defined protein-free medium. E. Aldermann et al., Proc.Natl.Acad.Sci.USA, vol 82, pp5771-5775, september 1985
HT-29 cells, protein free medium
62. 用于在存在饲养条件及无饲养条件下大规模培养人胚胎干细胞的永生饲养 (细胞)
Immortalized feeders for the scale-up of human embryonic stem cells in feeder and feeder-free conditions. Andre Choo, Jayanthi Padmanabhan, Angela Chin, Wey Jia Fong, Steve K.W. Oh, bioprocessing Technology Institute, Biommedical Sciences Institures, Singapore
Human embryonic stem cells, Immortalized feeders, Undifferentiated
63. 用琼脂糖凝胶电泳系统分离畜牧业兽用抗体
Agarose gel electrophoresis system for the separation of antibiotics used in animal agriculture. M.J. Salvatore et al., The Analyst 1993, vol. 118, no. 3, pp. 281-287
64. 从大鼠胶质瘤细胞株纯化糖蛋白血管内皮细胞线粒体
Purification of a Glycoprotein Vascular Endothelial Cell Mitogen From a Rat Glioma-Derived Cell Line. G. Conn et al., PNAS 1990, vol. 87, no. 4, pp. 1323-1327
65. 人骨髓干/祖细胞促进小鼠海马内源性神经干细胞的神经发生
Human stem/progenitor cells from bone marrow promote neurogenesis of endogenous neural stem cells in the hippocampus of mice. J.R. Munoz et al., PNAS 2005, vol. 102, no. 50, pp. 18171-18176

Culturing HEL 299 Cell Line on a Nunclon Δ Cell Culture Treated Surface

Introduction

NUNC Brand Nunclon. cell culture products are tested for cell growth and plating efficiency using several different cell lines. Nunclon. products are tested with two cell lines L929, HEL 299 or F2002 and one Primary Chick Embryo cellculture for monolayer formation, plus cell line V79 (-4) for cloning efficiency. The HEL 299 cell line is derived from embryonic lung tissue of a human male. It is a diploid fibroblast-like cell line initially developed for use in vaccine development.

This TechNote describes a procedure for culturing HEL 299 cell line on a Nunclon. treated surface.

Materials and Methods

- HEL 299 Cells (ATCC CCL 137)
- Minimum Essential Medium Eagle (MEM)
- Fetal Bovine Serum (FBS)
- L-Glutamine, 200 mM
- Sodium Pyruvate, 100 mM
- Lactalbumin Hydrolysate, 10% in Earle's Balanced Salt Solution
- Dulbecco's Phosphate Buffered Saline, 1X (without Ca^{2+} or Mg^{2+})
- EDTA, 0.02% Solution
- Trypsin Solution, 1X
- Antibiotic/Antimycotic Solution, 100X
- Crystal Violet or Methyl Violet, 0.1-0.4% in aqueous alcohol solution
- Reference lots

Prepare growth medium for HEL 299 cell line as follows:

MEM 1X	500.0 mL
FBS	58.0 mL
L-glutamine	5.8 mL
10% Lactalbumin Hydrolysate	5.8 mL
Sodium pyruvate	5.8 mL
Antibiotic/Antimycotic	5.8 mL
Total	581.2 mL

Culturing Procedure

1. Place culture vessels (samples and reference lots, as appropriate) in a laminar flow hood along with the Culture medium components which have been pre-warmed to 37°C.
2. Prior to harvesting, cells must be at least 75% confluent with good morphology. Aspirate medium and wash cells twice with 1X PBS before addition of EDTA solution.
3. Add EDTA solution to cover the Growth area completely
4. After EDTA solution is decanted, add trypsin to disaggregate the cells. Incubate culture vessels at 25°C or 37°C and monitor cell detachment under the microscope. Detachment time will vary.
5. After cells detach, add medium to stop trypsinization and to disperse the cells.
6. Transfer cells to a sterile conical tube and place on ice.

7. Determine cell quantity, e.g. Trypan Blue dye exclusion assay
8. Determine the number of cells required for each product to be Tested by multiplying the plating density by the surface area. Plating density for HEL 299 cell line is $2.0 \times 10^4/\text{cm}^2$.
9. Dilute cells into growth medium and Seed cell culture product.
10. Incubate cells in a 37°C incubator with 5% CO₂ for seven days to 2 form a confluent monolayer.
11. Decant medium. Add reagent alcohol, 95%, for 5 to 10 minutes for fixation, then decant.

Add crystal violet or methyl violet stain, 0.1-0.4%, to cover the surface for 5 to 10 minutes, then decant and wash with water.
12. Evaluate the monolayer when dry. (See Figure 1.)


Figure 1:

HEL 299 cells after seven days incubation at 37°C, cultured on a Nunclon Δ polystyrene surface, stained with 0.4% crystal violet. Calibration bar is 40 μ m.

Certification Results

When used for Nunclon. Certification, HEL 299 cell line results are evaluated as a percentage of surface coverage, per test sample.

- The average percent value must be within 10% of the values of the control products tested with HEL 299 cells.
- Cell growth must be consistent over the entire growth surface.

If these two conditions are met, product passes HEL 299 testing.

www.nuncbrand.com

Part of Thermo Fisher Scientific

Culturing L929 Cell Line on a Nunclon Δ Cell Culture Treated Surface

Introduction

NUNC™ Brand Nunclon Δ cell culture products are tested for cell growth and plating efficiency using several different cell lines. NunclonΔ products are Tested with two cell lines L929, HEL 299 or F2002 and one Primary Chick Embryo cell culture for monolayer Formation, plus cell line V79(-4) for cloning efficiency. L929 is a fibroblast-like cell line cloned from strain L. The parent strain was derived from normal subartaneous areolar and adipose tissue of a male C3H/An Mouse.

This TechNote describes a procedure for culturing L929 cell line on a Nunclon Δ treated surface.

Materials and Methods

- L929 Cells (ATCC CCL 1)
- Minimum Essential Medium Eagle (MEM)
- Bovine Calf Serum (BCS), iron supplemented or Newborn Calf Serum (NCS)
- L-Glutamine, 200 mM
- Non-essential Amino Acids (NEAA), 100X
- Dulbecco's Phosphate Buffered Saline, 1X (without Ca^{2+} or Mg^{2+})
- Trypsin Solution, 1X
- Antibiotic/Antimycotic Solution, 100X
- Crystal Violet or Methyl Violet, 0.1-0.4 % in aqueous alcohol solution
- Reference lots

Prepare growth medium for L929 cell line as follows:

MEM 1X	500.0 mL
BCS	57.0 mL
L-glutamine	5.7 mL
NEAA	5.7 mL
Antibiotic/Antimycotic	5.7 mL
Total	574.1 mL

Culturing Procedure

1. Place culture vessels (samples and reference lots, as appropriate) in a laminar flow hood along with the culture medium components which have been pre-warmed to 37°C.
2. Prior to harvesting, cells must be at least 75% confluent with good morphology. Aspirate medium and wash cells twice with 1X PBS before Trypsinization.
3. Add trypsin to disaggregate the cells. Incubate culture vessels at 25°C or 37°C and monitor cell detachment Under the microscope. Detachment time will vary.
4. After cells detach, add medium to stop trypsinization and to disperse the cells.
5. Transfer cells to a sterile conical tube and place on ice.
6. Determine cell quantity, e.g. Trypan Blue dye exclusion assay.
7. Determine the number of cells required for each product by Multiplying the plating density by the surface area. Plating density for L929 cell line is $1.5 \times 10^4/\text{cm}^2$.

8. Dilute cells into growth medium and seed cell culture product.
9. Incubate cells in a 37°C incubator with 5% CO for four days. A confluent monolayer of L929 cells should be formed.
10. Decant medium. Add reagent alcohol, 95%, for 5 to 10 minutes for fixation, then decant. Add crystal violet or methyl violet stain, 0.1-0.4%, to cover the surface for 5 to 10 minutes, then decant and wash with water.
11. Evaluate the monolayer when dry. (See Figure 1.)


Figure 1:

L929 cells after four days incubation at 37°C, cultured on a Nunclon Δ polystyrene surface, stained with 0.4% crystal violet. Calibration bar is 40μm.

Certification Results

When used for Nunclon Δ Certification, L929 cell line results are evaluated as a percentage of surface coverage, per test sample.

- The average percent value must be within 10% of the values of the control products tested with L929 cells.
- Cell growth should be consistent over The entire growth surface.

If these two conditions are met, product passes L929 testing.

www.nuncbrand.com

Part of Thermo Fisher Scientific

Culturing Primary Chick Embryo Cells on a Nunclon Δ Cell Culture Treated Surface

Introduction

NUNCTM Brand NunclonΔ cell culture products are tested for cell growth and plating efficiency using several different cell lines. NunclonΔ products are tested with two cell lines L929, HEL299 or F2002 and one Primary Chick Embryo (PCE) cell culture formonolayer Formation, plus cell line V79(-4) for cloning efficiency. Primary Chick Embryo cells are used to assess primary cell growth.

This TechNote describes a procedure for culturing cells isolated from chick embryos on a NunclonΔ treated Surface.

Materials and Methods

- Fertile Eggs (10 to 12 day gestation)
- Minimum Essential Medium Eagle (MEM)
- Bovine Calf Serum (BCS), iron supplemented or Newborn Calf Serum (NCS)
- L-Glutamine, 200 mM
- Non-essential Amino Acids (NEAA), 100X
- Dulbecco's Phosphate Buffered Saline, 1X (without Ca^{2+} or Mg^{2+})
- Trypsin Solution, 1X
- Antibiotic/Antimycotic Solution, 100X
- Crystal Violet or Methyl Violet, 0.1-0.4 % in aqueous alcohol solution
- Reference lot

Prepare growth medium for primary chick embryo cells as follows:

MEM 1X	500.0 mL
BCS	57.0 mL
L-glutamine	5.7 mL
NEAA	5.7 mL
Antibiotic/Antimycotic	5.7 mL
Total	574.1 mL

Culturing Procedure

A. Preparation of Chick Embryos

1. Place eggs, smaller (pointed) end up, in a sterile beaker.
2. Sterilize shell with 70% ethanol.
3. Pierce egg shell with sterile forceps and continue in a circular pattern to enlarge the opening until an intact circle of shell can be lifted from the Egg.
4. Extricate the embryo and decapitate. Place the body into a sterile 150x15 mm Petri dish.
5. Remove the limbs of the embryo with sterile scissors.
6. Place the remaining embryo in a 250 mL beaker and rinse twice with 1X PBS before mincing with sterile scissors.
7. Add 10 mL of trypsin (37°C, pH 7.5 to 8.0).
8. Stir the embryo-trypsin mixture on a magnetic stirrer for 45 minutes to 1 hour at 25°C.
9. Filter the mixture using sterile gauze Premoistened with 10 mL of 1X PBS.
10. Place the filtrate in a 50 mL conical tube and centrifuge for 10 minutes at 580 xg at 20°C.
11. After centrifugation, decant supernatant and resuspend cells in 10 mL of medium.

B. Culturing Procedure

1. Place culture vessels (samples and controls, as appropriate) in a laminar flow hood along with the culture Medium components which have been pre-warmed to 37°C.
2. Determine cell quantity, e.g. Trypan Blue dye exclusion assay. Each embryo yields approximately 3×10^8 cells.
3. Determine the number of cells Required for each product under test by multiplying the plating density by the surface area. Plating density for PCE cells is $1.5 \times 10^5/\text{cm}^2$.
4. Dilute the appropriate number of cells in growth medium and seed cell culture product.
5. Incubate cells in a 37°C incubator with 5% CO_2 for three days to form a confluent monolayer.
6. Decant medium. Add reagent alcohol, 95%, for 5 to 10 minutes for fixation, then decant. Add crystal violet or methyl violet stain 0.1-0.4%, to cover the surface for 5 to 10 minutes, then decant and wash with water.
7. Evaluate the monolayer when dry. See Figure 1.


Figure 1:

Primary chick embryo cells after three days incubation at 37°C, cultured on a NunclonΔ polystyrene surface, stained with 0.4% crystal violet. Calibration bar is 40μm.

Certification Results

When used for NunclonΔ Certification, primary cell growth results are evaluated as a percentage of surface coverage per test sample.

- The average percent value must be within 10% of the values of the control products tested with PCE cells.
- Cell growth must be consistent over the entire growth surface.

If these two conditions are met, product passes PCE testing.

www.nuncbrand.com

Part of Thermo Fisher Scientific

Culturing V79-4 Cell Line

on a Nunclon Δ Cell Culture Treated Surface

Introduction

NUNC™ Brand Nunclon Δ cell culture products are tested for cell growth and plating efficiency using several different cell lines. Nunclon Δ products are tested with two cell lines L929, HEL 299 or F2002 and one Primary Chick Embryo cell culture for monolayer formation, plus cell line V79-4 For cloning efficiency. V79-4 is a fibroblast-like cell line derived from the lung tissue of a male Chinese hamster. It has a relatively high plating efficiency and short generation time.

This Tech Note describes a procedure for culturing V79-4 cell line on a Nunclon Δ treated surface.

Materials and Methods

- V79-4 Cells (ATCC CCL 93)
- Minimum Essential Medium Eagle (MEM)
- Bovine Calf Serum (BCS), iron supplemented
- Fetal Bovine Serum (FBS)
- L-Glutamine, 200 mM
- Non-essential Amino Acids (NEAA), 100X
- Dulbecco's Phosphate Buffered Saline, 1X (without Ca^{2+} or Mg^{2+})
- Trypsin Solution, 1X
- Antibiotic/Antimycotic Solution, 100X
- Crystal Violet or Methyl Violet, 0.1-0.4% in aqueous alcohol solution
- Product Controls

Prepare growth medium for primary chick embryo cells as follows:

MEM 1X	500.0 mL
BCS	57.0 mL
L-glutamine	5.7 mL
NEAA	5.7 mL
Antibiotic/Antimycotic	5.7 mL
Total	574.1 mL

Culturing Procedure

1. Place culture vessels (samples and controls, as appropriate) in a laminar flow hood along with the culture medium components which have been pre-warmed to 37°C.
2. Prior to harvesting, cells must be at least 75% confluent with good morphology. Aspirate media and wash cells twice with 1X PBS before Trypsinization.
3. Add an appropriate volume of Trypsin to disaggregate the cells. Incubate culture vessels at 25°C or 37°C and monitor cell detachment under the microscope. Detachment time will vary.
4. After cells detach, add media to stop trypsinization and to disperse the Cells.
5. Transfer cells to a sterile conical tube and place on ice.
6. Determine cell quantity by the Trypan Blue dye exclusion assay.
7. Determine the number of cells required for each product by multiplying the plating density by the surface area. Plating density for V79-4 cell line is $6.0 \times 10/\text{cm}^2$.

8. Dilute the appropriate number of cells into growth media and seed cell culture product.
9. Incubate cells in a 37°C incubator with 5% CO_2 for six days to form distinct colonies.
10. Decant media. Add reagent alcohol, 95%, for 5 to 10 minutes for fixation, then decant. Add crystal violet or methyl violet stain, 0.1%-0.4%, to cover the surface for 5 to 10 minutes, then decant. And wash with water.
11. Evaluate cloning efficiency when dry. (See Figure 1).


Figure 1:

V79-4 cells after six days incubation at 37°C, cultured on a Nunclon Δ polystyrene surface, stained with 0.4% crystal violet. Calibration bar is 40 μm .

Certification Results

When used for Nunclon Δ Certification, V79-4 cell line results are evaluated as number of colonies per test sample.

- The average number of colonies must be within 15% or less of the average colonies observed in the control products tested with V79-4 cells.

If the previous condition is met, product passes V79-4 testing.

www.nuncbrand.com

Part of Thermo Fisher Scientific

Serum-free Transient Expression of Antibodies

Using 40 Tray Cell Factories

Peter Schlenke, Katrin Boetzel, Thilo Grob, Frank Kohne, Hans Hofer, Ulrike Brandle, Kurt Russ and Ralf Fehrenbach

Abstract

A serum-free transient expression method for antibodies in HEK Cells and a standardized purification protocol were developed. The methods described enable the fast and cost-effective production of hundreds of milligrams of different antibodies that can be used for test purposes. The N-Glycan profile of the antibodies obtained is comparable to that of antibodies expressed in stable CHO cell lines, confirming the suitability of the system to produce material for target validation studies.

Introduction

For target validation studies, significant amounts of antibodies are required. Because most of the biopharmaceuticals require specific post-translational modifications for proper bioactivity, e.g. Glycosylation, mammalian cells are mainly used as expression systems.

The generation of stable cell lines is expensive and time consuming and therefore restricted to a small number of antibody candidates, consequently transient expression of the antibodies offers a suitable method.

We have developed a transient serumfree expression method in HEK 293 cells using calcium-phosphate precipitation. The transfection of the cells is performed during cultivation in 40 Tray Cell Factories (CF 40), each with a culture area of about 25.000 cm². Due to the adherence of the cells, repeated harvesting is easily possible and results in a significant increase in the final product yield (see Fig.1). After a 3-step purification process a highly purified material is obtained. Analysis of the NGlycan structures released from the purified antibodies confirmed the suitability of the system to produce antibodies with characteristic glycosylation pattern.


Fig. 1 Harvesting of the antibody containing supernatant and replacement with fresh production medium using an Automatic Cell Factory

Results Transient expression

Two IgGs were expressed in HEK293 cells cultivated in 4 CF40. In each case 4 harvests were performed repeatedly after 3-4 days of cultivation. Production volume was 4 L/CF40, resulting in a total harvest volume of 64 L per IgG. The expression level reached up to 12 mg/L in serum-free culture conditions. For the two IgGs expressed a final product yield of about 530 resp. 430 mg was achieved. The results obtained are summarised in Tab. 1.

Tab. 1 Summary of the results obtained after large scale transient expression of two human IgGs

Harvest	1 st IgG			2 nd IgG		
	Titer [mg/L]	Volume [L]	Yield [mg]	Titer [mg/L]	Volume [L]	Yield [mg]
1	12	16	192	10	16	160
2	9	16	144	9	16	144
3	8	16	128	6	16	96
4	4	16	64	2	16	32

□ 528 □ 432

Purification of the IgGs

The IgGs were purified from the 64 L cell culture supernatant using a standardised 3-step protocol after concentration by ultrafiltration. SDS PAGE analysis, isoelectric focusing and HPLC analysis (data not shown) confirmed, that the purification strategy resulted in a highly purified product (see Fig. 2).


Fig. 2 SDS PAGE analysis and isoelectric focusing. The IgG was analysed during and at the end of the purification process. After separation the IgG was visualized by silver staining.

Analysis of the N-Glycan structures

The main N-Glycans detected after desialylation were complex-type biantennary structures without terminal Gal (36.6%), with 1 terminal Gal (44.9%) and with 2 terminal Gal (15.1%). Only about 2% of the NGlycans were sialylated. The analysis confirmed the ability of the transiently transfected HEK293 cells to synthesise N-Glycan structures similar to those obtained after stable expression in e.g. CHO cells.


Fig. 3 HPAEC PAD profile of the desialylated NGlycans released from IgGs after transient expression in HEK293 cells. The structures detected are summarised in Tab. 2.

Tab. 2 Desialylated N-Glycan structures found on an IgG after transient expression in HEK293 cells.

Peak	Retention [min]	Structure Code	Area [%]	Structure
1	11.23	CoreF1GN1G050	0,9	
2	14.43	CoreF1GN2G050	36,6	
3	17.60	CoreF1GN2G130	44,9	
4	20.37	CoreF1GN2G250	15,1	
5	23.13	not determined	1,35	
6	24.30	not determined	1,18	

GalNAc: ■ Man: ● Gal: ■ Fuc: ◆

Conclusion

- By transient expression of two human IgGs in 4x40 Tray Cell Factories, a final product titer of 530 respectively 430 mg was obtained.
- A 3-step purification process using a standardised protocol resulted in a highly purified product.
- The N-glycosylation pattern of the IgGs, after transient expression in HEK293 cells, were similar to IgGs synthesised in stably transfected CHO cell lines.

Application Note


Production of Poliovirus using Vero Cells Cultured in Nunc Cell Factories, T-Flasks and Roller Bottles with HyQ SFM4MegaVir

Overview

Biopharmaceutical markets demand that many products be made on a large scale. As researchers develop new products, they are required to take a process from bench scale to production scale. Those who have gone through the scale up process realize that simply putting more of the same ingredients in a larger container may not produce the same results. Also, container scale up may not be feasible. For instance, if a researcher required a 100-fold increase in virus, the surface area required for a T-flask would need to be increased from 75 cm² to 7,500 cm². Imagine the size of a T-7500 flask (7,500 cm² or 8 square feet!).

Many approaches exist for cell culture scale up and the production of virus. Three cell culture vessels (Nunc Cell Factories, Corning T-75 flasks, and Corning roller bottles) are compared here in terms of why and how each is used, protocol differences, and virus production performance.

The biopharmaceutical industry is moving away from the use of animal derived products used in the manufacturing of products. The MegaVir medium and non-mammalian trypsin replacement used in this study were selected to follow this trend. MegaVir is a protein-free medium that is used without the addition of animal serum. MegaVir has been shown to sustain anchorage-dependent VERO cell cultures over long periods. Vaccine manufacturers use VERO cultures to produce viruses such as the poliovirus used in this study.

Materials and Methods

VERO cells (ATCC CCL-81, African Green Monkey, adult kidney, epithelial) were adapted from serum-containing conditions to MegaVir medium (HyClone) and carried over 20 passages. Cells were then cultured in Corning T-75 cell culture flasks and scaled up to Nunc Cell Factory 10-layer unit (CF10) and 850 cm² roller bottles (Corning) on a Belco Technology Cell-Production Roller Apparatus for virus propagation studies. Surface area and medium volumes for each vessel type are outlined in table 2.

Cells were cultured at an initial density of 8,000 cells/cm², cultured at 37°C in a 5% CO₂ / 95% air environment, and observed daily under a microscope until at least 2 90% confluency was achieved. Medium was then removed from the cell culture vessel, cells were rinsed with phosphate buffered saline solution (PBS), and nonmammalian trypsin replacement was applied at 6-10 mL per 75-cm² area. Cells were observed until all cells had detached from the culture surface (typically 15-20 minutes), then cell counts and viability were performed, recorded, and analyzed for doubling times. The non-mammalian trypsin replacement was not removed from the culture. New vessels were seeded for scale up at 8,000 cells/cm².

Upon reaching 90% confluency the cultures were inoculated with poliovirus at an approximate MOI of 0.1. The culture was then incubated at 37°C for 72 h, followed by aseptic removal and storage of the supernatant. Any cells remaining attached to the growth surface of the CF were removed using non-mammalian trypsin replacement, and added to the supernatant. In order to lyse all cells the supernatant was completely frozen at -70°C then thawed completely at 37°C. This process was repeated twice to lyse cells and release virus particles. Lysate was then collected and centrifuged to pellet and remove cellular debris. The supernatant containing virus was used to prepare 10 log dilutions of each sample. Titrations for each 10 clarified sample were performed in 96 well micro-titer plates using TCID₅₀ assays 50 described previously.

Specific protocols and advantages for each vessel are outlined in Table 1.

Table 1	Market Application	Additional Equipment Required	Gas Transfer Mechanism	Fluid Movement
T-Flask	Research	Incubator/ Pipetting System	Cap	Pipette
Roller Bottles	Research and Industry	Incubator, Pipetting System, Roller Apparatus	Cap	Pipette
Cell Factories	Research and Industry	Incubator/ Pipetting System	Filter	Pipette and Gravity Feed via Tubing

Results and Discussion

Vero cells in this study reached confluency normally in all vessels tested with good surface attachment, normal morphology, and even distribution of cells. All doubling times and peak cell densities at the times of passage fell in a normal range for protein-free cell culture conditions. Constant doubling times (approximately 40 hours) and good viability were maintained.

No detrimental affects were observed from the use of non-mammalian trypsin replacement. Handling was reduced since there is no need to centrifuge cells to remove the cell detachment solution before reseeding a culture vessel.

Scaling from a T-75 to a CF1 to a CF10 unit was possible in two passages. The 85-fold increase of surface area/ cell mass is similar to scaling up suspension cultures from a 250 mL volume shaker flask to a 15-L bioreactor in two passages.

Virus Production

T-flasks are typically used for research applications, developing cell lines, etc., but not convenient for large scale applications since considerable time is required to handle the large number of vessels required. Roller bottles are an improvement over T-flasks since a greater surface area is available compared to T-flasks while using less media. A disadvantage of the roller bottle method is the requirement for a roller bottle apparatus to rotate the cylinders, supplying nutrients and gas exchange to the entire surface. Also, the apparatus greatly increases the footprint or incubator surface area required for operation. Physical space is often a highly coveted property in labs and production facilities.

A great advantage of the Nunc Cell Factory is its small footprint and its scalability from R&D through large-scale production (see table 2). A short learning curve occurs as operators learn to dispense media, trypsin, and PBS into and out of the vessel (see manufacturer's instructions for a detailed description). Larger Cell Factories (CF10 and CF40) would require the use of a non-mammalian trypsin replacement. Larger vessels mean larger volumes of media that make centrifugation necessary to remove cell-damaging trypsin. It is possible to leave a non-mammalian trypsin replacement in cell passages. This method decreases the labor involved as well since centrifugation is not required.


Table 2. Comparing surface area, reagent volume, total footprint, labor, and additional equipment requirements for T-flasks, roller bottles, and Nunc Cell Factories.

Criteria	Cf10	T-75	Roller Bottle
Surface Area (cm ²)	6,320	75	850
Media Volume/Vessel (mL)	2,000	50	200
Units Required for 6,320 cm ² area	1	85	8
Total Media Required for 6,320 cm ² area (mL)	2,000	4,250	1,600
Total Footprint (cm ²)	670	2,670*	4,250**
Total Non-Mammalian Trypsin Replacement (mL)	250	510***	160**
Total Volume PBS (mL)	250	850***	160****
Labor Hours to Harvest	0.5	2.0	0.5
Virus Titer (TCID ₅₀ / mL)	1.26 × 10 ⁸	2.61 × 10 ⁸	1.13 × 10 ⁹
Total Virus Titer (TCID ₅₀)	2.52 × 10 ¹¹	1.31 × 10 ¹⁰	2.26 × 10 ¹¹

* Based on stacking T-75 flasks four high
** Based on Bellco Technology Cell Production Roller Apparatus
*** Based on 6 mL non-mammalian trypsin replacement and 10 mL PBS/flask
**** Based on 20 mL non-mammalian trypsin replacement and 20 mL PBS/flask

Summary

Protein-free cell culture is becoming more commonplace in biopharmaceutical production. The combination of HyQ SFM4MegaVir and the non-mammalian trypsin replacement suit this role and perform well in cell growth and virus production. The medium contains no animal derived components and the monolayer-detaching enzyme used is a non-trypsin, non-mammalian reagent that is well suited for biopharmaceutical production methods.

Since HyQ SFM4MegaVir performs equally well in all systems tested, the choice of which vessel to use will depend on interests in time, convenience, and space savings. Roller bottles may outperform Cell Factories on virus production per mL of medium used, but footprint and labor savings may outweigh the media savings benefit.

Ordering Information

Vendor	Description	Part Number	Volume	Delivery System
HyClone	HyQ SFM4MegaVir	Sh30552	500mL-1000mL	PETE Bottles
			5L-500L	CX5-14 BPC
Vendor	Description	Part Number	Number of Trays	Units per Case
Nalge Nunc	CF-10 Nunc Cell Factory	164327	10	2
		170009	10	6

Scaleable, controlled growth of adherent cells in a disposable, multilayer format

Extending Cell Factories with uniform gas-distribution

Edwin Schwander, Hans Rasmusen, NUNC A/S

History

Until the 1970's, the larger scale production of adherent cells was mainly performed in bottles, where the distribution of cells and medium was achieved by a rotating motion. This technology has the advantage of being simple and easy in operation, achieving an increased volume-surface ratio of about $0.3 \text{ cm}^2/\text{cm}^3$ (standard Roller Bottle). Bottles contain a sufficient volume of gas with reasonable gas transfer rates, but are limited in terms of size, to about 4000 cm^2 per unit. They are cumbersome to operate manually, when scaling up is needed.


In 1975, Nunc A/S began to produce a new generation of devices, in collaboration with Bioferon (Renschler Biotechnology). This device, called Cell Factory, consists of trays with a surface of 632 cm^2 each. Stacked trays of up to 40 layers (with a total surface of $25,280 \text{ cm}^2$ per unit) achieve a surface-volume ratio of $0.5 \text{ cm}^2/\text{cm}^3$. This creates a format, for using a flask-like system from research applications up to industrial production scale.

In order to industrialize the use of such units, a robotic system called ACFM (1) (Automatic Cell Factory Manipulator) was developed to physically manipulate 4 interconnected 40-layer units simultaneously.

In this way a surface of $101,120 \text{ cm}^2$, representing the equivalent of about 120 standard Roller Bottles, could be handled as one closed unit. With limited space requirements (less than 150 m^3), batches of several thousand liters can be produced with little labour involved. Gradually, a whole line of equipment, from handling/supports for single units (2+3) to shakers (4), rack-systems (5), incubators (6), warm-room-shelves (7), was developed.

Nunc also has a series of new products under development such as simple gassing stations (8), microscopes (9) and motion-generators (10).

This system found its way into the pharmaceutical community for the production of vaccines and proteins (polio-, rabies-, Hep-A-vaccine, EPO, interferons, hGH and others products) at intermediate to very large scale production of several thousand liter per batch.

Key aspects of multilayer format


The Cell Factory as a device compared to Roller Bottles offers the significant advantages that with a 120 x surface, compared to standard Roller Bottles, labour time and contamination risks are drastically reduced.

This compact format is realized mainly by the reduction of the air volume above the culture.

In recent years, the use of different, partly modified, cell lines created a situation, where the availability of oxygen became a limiting factor in some cultures.


Gassing the devices, by just putting airfilters on both channels, is possible, but delivery of the gasses is uneven within the device due to negative back pressure (see figure 1).

figure 1

In order to overcome this problem, Nunc created a patented product extension, where a back pressure tube is inserted in one of the filling channels (see figure 2). This tube with 0.4 mm pinch holes for each layer, creates sufficient back pressure to insure that the gasses are distributed equally over the individual layers at a moderate flow rate of 50 mL per tray.

figure 2


In an experimental set-up, oxygensensors were placed in tray 1, 4, 7 and 10 (counted from the top) of a gassed 10- layer Cell Factory and a 10-layer unit of a competitor. Units were filled with 2 liter of a buffer and oxygen was depleted by extensive flushing with nitrogen.

Then units were gassed with a constant flow (0,5l/min) of atmospheric air and oxygen saturation was measured over a period of 90 minutes.

In figure 3 the results are shown. While with the gassed unit an almost identical pattern of saturation with oxygen can be seen, the competitor product shows a significant difference in saturation in time, indicating clearly a non-identical cell culture environment.

figure 3

First results comparing active gassed to non-gassed cultures

Nadia De Bernardi, Luisa Nolli, Areta International srl

In order to investigate the effects of active gassing, the following studies were initiated at Areta International srl (Gerenzano, Italy):

Cells were cultivated under active gassed and non-gassed conditions, monitoring glucose/lactate and comparing cell density.

First results obtained from these experiments indicate the following:

- Higher growth rate was observed when gassing CF's
- Cell densities reached were on average 20-40% higher under gassing
- Expression of proteins (IgG) was higher under active gassing conditions
- Faster adaptation to serum free medium (SFM) was observed in experiments

In each experiment, different cell lines, both suspension and anchorage dependent, were cultivated. In all experiments 1×10^4 cells/cm² in 400 mL of medium were used for seeding the devices.

The cells were maintained in culture for 4 days. During the culture period, the medium was sampled daily from each CF to test viability, cell number and at the end of the culture time, the productivity.

In gassed CF's a pump was set to pass 2 litres of gas into the CF's for one minute every four minutes. The non gassed CF's were maintained in an incubator with 5% CO₂.

figure 4


figure 5



figure 6

The figures below show data and results obtained.

In summary gas is beneficial for:

- Cell expansion (figure 4: cell density is increased by up to 73% when seeded in Cell Factories with gassing)
- Productivity (the production is 58% higher in a gassed Cell Factory) (figure 5)
- Serum free adaptation: the cell growth is faster in CF with gas in SFM. (figure 6)

Many disposable systems for cultivating adherent cells have appeared on the market during the last 30 years, but few of them reached a broad public, and some did not progress beyond the proto-type stage. Nunc's Cell Factory has proven to be a viable format for the production of adherent cells. The fact that competitors recently started to copy the product, both confirms the security of the market position, and shows that the Cell Factory format has reached a mature status.

The key factor when deciding on a system is an assurance of the commitment of the supplier to support its product over a long period of time. In this respect, over 3 decades Nunc has proven to be strongly committed to Cell Factory products.

Continuous ongoing improvements and customized product extensions are part of our philosophy.

Acknowledgements

The authors want to thank Vibeke Rowell, Ron Main and Birgitte Raldorf for reading and revising the manuscript.

Thermo Scientific Nunc Cell Factory System extractable review

How product design influences an extractable profile

Stephanie Carter, Application Scientist; Joseph Granchelli, PhD, Manager, Applications and Technical Support Thermo Fisher Scientific, Rochester, New York

Key Words

Cell Factory System, extractables, UV-cured adhesive, ultra-sonic welding, product design, risk, safety

Introduction

Biopharmaceutical manufacturers have an obligation to ensure the safety and efficacy of the products they produce. This requirement drives the selection of products and materials used in manufacturing operations and define the necessary risk assessment of those materials.

One of the areas of concern are the chemicals that can be extracted from cell culture vessels and other containers used in the storage and manufacture of the pharmaceutical. To mitigate these risks, biopharmaceutical manufacturers often use extractables data to perform a risk assessment to determine if they need to take action to reduce risk or if further studies are required. These studies are performed using conditions which are exaggerated to ensure that the list of potential extractables is comprehensive.

Suppliers can assist manufacturers by incorporating an understanding of extractables into the design of their products. Not only is the selection of materials important, but designing products in a way that minimizes the number of materials can reduce the amount of time it takes to assess

risk and help save manufacturers' time. The more individual materials you incorporate into a product design, the greater the potential to increase the number of extractables observed.

Thermo Scientific™ Nunc™ Cell Factory™ systems reflect an understanding of the impact of product design on the potential to generate extractables. The Nunc Cell Factory System minimizes the number of materials that could generate extractables by relying on ultrasonic welding to fuse individual layers of the Cell Factory system together. A Competitor uses an adhesive to bond individual layers Thermo Scientific Nunc Cell Factory System extractable review - How product design influences an extractable profile Stephanie Carter, Application Scientist; Joseph Granchelli, PhD, Manager, Applications and Technical Support Thermo Fisher Scientific, Rochester, New York of their multilayer vessel, creating the potential for materials from that adhesive to migrate from the vessel into the cell culture media. The Nunc Cell Factory system has no adhesives; therefore, extractables associated with UV-cured adhesive are absent from this system.

An extractable study was conducted to demonstrate the potential for an adhesive to leach from the Competitor's multilayer cell culture vessel in comparison to the adhesive-free Nunc Cell Factory and Thermo Scientific™ EasyFill™ Cell Factory™ systems.

Materials

Nunc Cell Factory 4-Layer Standard System;

Catalog number 140004; Lot number 1099245; Dates of Manufacture: July 03 to 09, 2013

Nunc EasyFill Cell Factory 4-Layer System; Catalog number 140360; Lot number 136652; Dates of Manufacture: November 2013

Competitor 5-layer vessel; Date of Manufacture: December 17, 2013

Methods

The Nunc Cell Factory Standard System, Nunc EasyFill Cell Factory System, and a Competitor's multilayer vessel were filled with 800mL of phosphate buffered saline (PBS, pH 7.4), or 20% isopropyl alcohol (IPA), and were incubated for 20 days at 40° C to exaggerate conditions of use (n=1 per test condition¹). PBS was chosen because it is very similar in ionic strength and polarity to cell culture media and to human biological fluids, and as a result should represent a physiologically relevant solution. IPA is slightly more aggressive than PBS in terms of its ability to dissolve certain chemicals.

All extractions and analysis were performed at a thirdparty laboratory.

The extractions were analyzed by Direct Injection GC/MS, Headspace GC/MS, and LC/MS for semi-volatile organic compounds, volatile organic compounds, and non-volatile organic compounds, respectively.

No assessment was made of the toxicological implications of the presence or quantity of the extractables detected in any

¹ We expect from reasonable scientific principles that we should find the signature of an adhesive in an extractables report for the competitive product and not in the Cell Factory product. The sample size of one is confirmatory.

Table 1. Extracted Materials Identified in the Nunc Cell Factory and EasyFill Cell Factory Systems

Nunc Cell Factory 4-Layer Standard System		Nunc EasyFill Cell Factory 4-Layer System	
IPA Extracts	PBS Extracts	IPA Extracts	PBS Extracts
Benzaldehyde	None Detected	Benzaldehyde	None Detected
Styrene		One Unknown Compound	
Two Unknown Compounds			

Table 2. Extracted Materials Identified in the Competitor 5-layer vessel

Competitor's Multilayer Vessel	
IPA Extracts	PBS Extracts
(1-hydroxycyclohexyl)phenyl-methanone*	(1-hydroxycyclohexyl)phenyl-methanone*
2,4,6-Trimethylbenzaldehyde	N,N-dimethyl-2-propenamide*
Benzaldehyde	Four Unknown Compounds
Exobornyl Acetate*	
Isobornyl Alcohol*	
N,N-dimethyl-2-propenamide*	
Silicone	
16 Unknown Compounds	

* Known components of adhesives

vessels. No suggestion of suitability or unsuitability for any application should be implied based on these analyses.

Results and Discussion

The extractables profile of the Nunc Cell Factory Systems differed as expected (Table 1). Additional extractables were detected in the competitive device. The compounds detected in the competitor's vessel are consistent with what would be expected from an adhesive (Table 2). No adhesive is used in the construction of the the Nunc Cell Factory systems. As expected, no extractables from adhesive were detected in the Cell Factory products.

Conclusions

Specific design features have trade-offs. In the present example, the use of an adhesive to bond layers of a multilayer tray together imparts certain advantages in manufacturing and makes these joints more resistant to breakage under certain conditions. However, the presence of this additional material also carries with it other consequences that may or may not introduce a risk to the products produced in that device. This study demonstrates that the presence of an additional material alters the extractable profile of that device by generating additional extractables. The magnitude and significance of the risk posed by these additional extractables must be assessed where that product is used to produce a biopharmaceutical. Risk assessment can be a timeconsuming process and can lead to the decision to perform expensive follow-up tests, which can delay the introduction of a product to the market. Minimizing the number of materials in the product design is one way suppliers can assist customers in simplifying the risk assessment process.

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