

ProteinSEQ™ Protein A Quantification Kit

Catalog Number 4469343

Pub. No. MAN0013471 Rev. A.0

Box Part No.	Contents	Item Part No.	Quantity	Volume	Storage conditions
A25001	Repligen Protein A Standard (1 mg/mL)	4469347	1	0.1 mL	-20°C
	Protein A 5' Probe	4469376	1	0.7 mL	
	Protein A 3' Probe	4469381	1	0.7 mL	
	2X Protein A Sample Prep Reagent	4485137	1	10 mL	
	ProteinSEQ™ Ligation and Assay Mix	100024250	1	0.6 mL	
	ProteinSEQ™ Ligase	A25020	1	55 µL	
A25002	Fast Master Mix, 2X	4448622	1	6 mL	4°C
	Wash Buffer	4469398	1	100 mL	
	Protein A Capture Beads	4469368	1	1.2 mL	
	ProteinSEQ™ Elution Buffer	100024368	1	7 mL	
	ProteinSEQ™ Diluent	4469405	1	41 mL	

Note: For safety and biohazard guidelines, refer to the “Safety” appendix in the *ProteinSEQ™ Protein A Quantification Kit User Guide* (Pub. no. MAN0013524). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

The information in this guide is subject to change without notice.

DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/techresources or email techsupport@lifetech.com

thermofisher.com/lifescience