

ProteinSEQ™ CHO HCP Quantification Kit

Catalog Numbers A27601

Pub. No. MAN0010565 Rev. B.0

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *ProteinSEQ™ CHO Host Cell Protein Quantification Kit User Guide* (Pub. No. MAN0010806). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product Description

Use the ProteinSEQ™ CHO HCP Quantification Kit to quantify the Chinese Hamster Ovary (CHO) Host Cell Proteins (HCPs) present in your bioprocess sample(s).

Contents

Table 1 ProteinSEQ™ CHO HCP Quantification Kit (Cat. No. A27601)

Contents	Amount	Storage
Box (Cat. No. A25500)		
CHO HCP Standard (15.6 µg/mL)	0.25 mL	–25°C to –15°C Store at 2–8°C after thawing. ^[1]
CHO HCP 5′ Probe	0.35 mL	–25°C to –15°C
CHO HCP 3′ Probe	0.35 mL	
ProteinSEQ™ Ligation and Assay Mix	0.6 mL	
ProteinSEQ™ Ligase	55 µL	
Box (Cat. No. A25499)		
Fast Master Mix, 2X	6 mL	2–8°C
Wash Buffer	100 mL	
CHO HCP Capture Beads	5.5 mL	
ProteinSEQ™ Elution Buffer	7 mL	
CHO HCP ProteinSEQ™ Diluent	50 mL	

^[1] After thawing, do not re-freeze. Store at 2–8°C for up to 1 month.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Manufacturer: Thermo Fisher Scientific | 7 Kingsland Grange | Warrington, Cheshire WA1 4SR | United Kingdom

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0010565

Revision	Date	Description
B.0	20 April 2018	Update template and legal information.
A.0	17 Nov 2014	New document.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.